

FENDERWEDGE®

Interproximal Tooth Shields

Pre-separate and Protect™



“Even those with steady hands and stellar technique occasionally and inadvertently nick the proximal tooth or restoration as they are preparing the adjacent tooth. I know I have! The Fender Wedge is a simple device that I wish I would have had the opportunity to use since the day I started practicing dentistry. It would have saved me much aggravation and time. It is now a “must-have” item during all my restorative procedures.”

David S. Hornbrook, DDS, FAACD
International Lecturer and Private Practice,
San Diego, CA

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations.

FenderWedge is a combination of a steel plate and a plastic wedge. inserted into the interproximal space, it protects the tissue and separates the teeth, thus simplifying the application of a matrix.

Basic Technique

1. Place FenderWedge prior to preparation to pre-separate and protect.



2. Use a FenderWedge large enough to both pre-separate the teeth and firmly hold the protective shield in place even after contact is broken.



3. Prepare the tooth with a conservative Class II preparation.

4. Remove the FenderWedge when the preparation is complete. Separation created by the FenderWedge will facilitate easier matrix band placement.

FenderWedge® Kits

Description

Complete Assorted Kit

Item **Qty.**
TSK01 80

Refills	Item	Standard		Bulk	
		Item	Qty.	Item	Qty.
X-Small Purple	TSPR-M	30	TSPR	60	
Small Orange	TSOR-M	30	TSOR	60	
Medium Green	TSGR-M	30	TSGR	60	
Large Yellow	TSYL-M	30	TSYL	60	




Large (yellow)
#TSYL


Medium (green)
#TSGR


Small (orange)
#TSOR


X-small (purple)
#TSPR

IMPORTANT NOTICE: FenderWedge is NOT a sectional matrix system nor does it replace the use of a matrix band. It protects the adjacent tooth during preparation and provides a pre-wedging force.



For more information go to www.garrisondental.com. Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

To order: Call toll free **888.437.0032**
or online at www.garrisondental.com



Mention Code
IN0216TS