# Composi-Tight® 3D XR

**User instructions** (ENGLISH)

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# Composi-Tight® 3D XR

## User instructions.

The Composi-Tight® 3D XR Sectional Matrix System is designed to create tight, anatomical contacts on Class II composite restorations. The contoured sectional matrix bands create proper tooth anatomy while two styles of rings create the required tooth separation and band adaptation for a tight, natural contact.

For Patent information see [www.garrisdental.com/patents](http://www.garrisdental.com/patents)

Non-sterile, RX only, always use a rubber dam.

## Important Cleaning and Sterilization Instructions:

**Composi-Tight 3D & 3D XR™ Soft Face™ rings** remove all visible debris after use and prior to sterilization. DO NOT use solutions containing acetone or ketone. Steam autoclave sterilization only.

**Sectional Matrix Bands:** Matrix bands are single use disposable items. Sterilization prior to use is recommended and can be performed using steam autoclave.

**Instruments & Forceps:** - Instruments may be dry heat sterilized, steam autoclaved or cold sterilized.

**Wedges:** Wedges are single use, disposable items. Spray disinfect only. Do not autoclave.

**For complete sterilization instructions go to:** [garrisdental.com/sterilization](http://garrisdental.com/sterilization)

## Technique for Class II restorations utilizing 3D or 3D XR Soft Face™ Rings

1. **Place rubber dam.** Isolation of the area helps to prevent contamination from blood and saliva. As with many dental procedures, aspiration of small parts is a valid concern.
2. **Place FenderWedge** prior to preparation to pre-separate and protect the adjacent teeth. (figure 1) Prepare tooth with conservative Class II preparation. Once the preparation is complete, remove the FenderWedge and place your matrix.
3. **Place contoured sectional matrix band.** Select a matrix band that most closely approximates the occlusogingival height of the tooth. Over-approximate tooth curvature buccolingually by rolling the contoured band with your fingers (Fusion Full Curve matrix bands do not require rolling to wrap around the tooth). This will help keep the matrix band ends from interfering with the Soft Face™ 3D/3D XR Ring placement. The band should be oriented with the convex edge towards the gingival margin of the tooth. (figure 2)
4. **Insert wedge.** For optimal contour, gingival seal and tooth separation, use a Garrison wedge. Keep a finger on the matrix band to prevent the wedge from dislodging the band during insertion. Wedge firmly, a high viscosity resin packed into the preparation may displace the band away from the tooth. (figure 3)
5. **Apply the Soft Face™ 3D/3D XR ring** retainer by spreading it with the Composi-Tight Ring Placement forceps and placing it over the wedge (figure 4). NOTE: For secure gripping of ring tips to surface, dry teeth thoroughly before placement.
6. **Very Important:** push the matrix band firmly into contact with the adjacent tooth in the desired contact area.
7. **Use your preferred etching,** disinfecting and bonding technique following manufacturer's instructions. Apply your preferred composite resin following manufacturer's instructions.
8. **Remove the ring, wedge and band.** Because of the superior contact achieved with Composi-Tight 3D/3D XR, removal of the band may require a Composi-Tight Matrix Forceps (model #BFA-1) or hemostat. This is especially true with particular bonding agents that preferentially bond to metal. For easier removal choose Slick Bands™.
9. **Contour and polish the restoration.**

## Stacking of Rings & MOD Restorations

Both the 3D XR Soft-Face™ Ring (blue tip) and the 3D Soft-Face™ Ring (orange tip) are angled to allow stacking of rings for performing M.O.D.s and multiple-tooth restorations. Only one ring may be stacked over the other, and in any combination. (figure 5)

\*(Garrison part number AUMRDF-100 Ring Placement Forceps will hold the rings more securely and provide ample force to open the rings. NOTE: Most rubber dam clamp forceps will not open the 3D Rings far enough to allow placement.)

figure 1

figure 2

figure 3


figure 4

figure 5

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