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1698r0715

WedgeWands®

English: Important and Helpful Wedge Placement Tips:

Wedge Selection: Select a wedge size that will penetrate completely through the interproximal space. If only the tip is inserted, the wedge may shift or bend during wand removal.

Bend: IMPORTANT! Bend the Wedge slowly to a full 90° angle. This creates a stress in the plastic that is essential for twist-off of the wand, while improving access in hard-to-reach areas. A wedge bent only 45° will not separate from the wand easily.

Alternate Method: For greater insertion force, drive the wedge straight into the interproximal space without pre-bending; bend to a full 90° before attempting to twist off the wand.

Insert: Drive wedge firmly and completely through the interproximal area, gripping the wand as close to the wedge as possible for maximum leverage. If only the tip is inserted, the wedge may shift or bend when the wand is being twisted off. Full insertion will also provide maximum gingival adaptation and tooth separation.

Twist: After proper bending and insertion, a simple 360° twist will separate the wand from the wedge.

Alternate Method: The wand may remain attached to the wedge if it does not interfere with the work area. This allows easy removal of the wedge when the restoration is complete.

Insertion with traditional instruments:

Clinicians who prefer the unique shape and design of Wedge Wands, but prefer to insert with hemostats, may elect to remove the wedge from the wand prior to use.

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