






# Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials	Pros and Cons
<b>Back-to-Back</b> 	<b>BEST CHOICE:</b> Quad Ring, Quad Wedge, Firm Band	<b>Pros:</b> Ring driver tip and split wedge work together. Firm Bands prevent deformation. Reduced patient chair time. <b>Cons:</b> Learning curve.
<b>Wide Prep</b> 	<b>BEST CHOICE:</b> Strata-G Wide Prep Ring, Firm Band <b>OPTION 2:</b> 3D Fusion Wide Prep Ring, Firm Band <b>OPTION 3:</b> Reel Matrix	<b>Pros:</b> Smaller ring diameter. The only sectional ring that can perform on wide preps. <b>Cons:</b> Limited application. Not recommended for basic Class II. <b>Pros:</b> Similar performance to Strata-G. <b>Cons:</b> Larger ring diameter. <b>Pros:</b> For teeth too broken down for a sectional matrix system. <b>Cons:</b> Less ideal contour and contact compared to a sectional matrix solution.
<b>Deep Prep</b> 	<b>BEST CHOICE:</b> Quad Ring, Quad Wedge, and Firm Band <b>OPTION 2:</b> Strata-G Ring, Strata-G Wedge, and Firm Band	<b>Pros:</b> Split tip wedge sits deep in sulcus and is less likely to be pushed up by the gingiva. <b>Cons:</b> N/A <b>Pros:</b> Strata-G wedge sits deeper than most wedges. <b>Cons:</b> Strata-G wedge does not sit as deep as Quad.
<b>Class III</b> 	<b>BEST CHOICE:</b> Anterior Fusion <b>OPTION 2:</b> Vari-Strip	<b>Pros:</b> Natural emergence profile and interproximal contours. <b>Cons:</b> Learning curve. Shape may not be appropriate for all teeth. (i.e. bulbous). Not cure through. <b>Pros:</b> Simple. Cure through. <b>Cons:</b> Susceptible to deformation during use.
<b>Class IV</b> 	<b>BEST CHOICE:</b> Anterior Fusion <b>OPTION 2:</b> Vari-Strip	<b>Pros:</b> Natural emergence profile and interproximal contours. <b>Cons:</b> Learning curve. Shape may not be appropriate for all teeth. (i.e. bulbous). Not cure through. <b>Pros:</b> Simple. Cure through. <b>Cons:</b> Susceptible to deformation during use.
<b>Class V</b> 	<b>BEST CHOICE:</b> Blue View Cervical	<b>Pros:</b> Compresses composite into the restoration. Limited finishing required. <b>Cons:</b> N/A




# Matrix Selection Guide

## Special Cases

Procedure	Recommended Materials	Pros and Cons
<b>Deep Margin Elevation</b> 	<b>BEST CHOICE:</b> ReelMatrix, Deep Margin Elevation Band (P/N RM500H) <b>OPTION 2:</b> Deep Margin Elevation Tofflemire Band (P/N RM500H) with a standard tofflemire retainer	<b>Pros:</b> Ability to seal very deep lesions. Lightweight, retainerless. <b>Cons:</b> Will not restore whole tooth. Step one of a two step process. <b>Pros:</b> Ability to seal very deep lesions. <b>Cons:</b> Retainer can torque the band out of position.
<b>Last Tooth in Arch/ Missing Adjacent Tooth</b> 	<b>BEST CHOICE:</b> ReelMatrix <b>OPTION 2:</b> Tofflemire Slickbands	<b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. Retainerless bands allow for simultaneous multi-tooth restoration. <b>Cons:</b> N/A <b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. <b>Cons:</b> Retainer can torque the band out of position. Cannot perform simultaneous multi-tooth restorations.
<b>Core Buildup</b> 	<b>BEST CHOICE:</b> ReelMatrix <b>OPTION 2:</b> Tofflemire Slickbands	<b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. Retainerless bands allow for simultaneous multi-tooth restoration. <b>Cons:</b> N/A <b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. <b>Cons:</b> Retainer can torque the band out of position. Cannot perform simultaneous multi-tooth restorations.
<b>Pedo</b> 	<b>BEST CHOICE:</b> Strata-G Universal Ring <b>OPTION 2:</b> ReelMatrix	<b>Pros:</b> Retention capability on shorter teeth. <b>Cons:</b> N/A <b>Pros:</b> Speed. <b>Cons:</b> Less ideal contour and contact to a sectional matrix solution.
<b>Distal of Canine</b> 	<b>BEST CHOICE:</b> Strata-G Universal Ring <b>OPTION 2:</b> 3D Fusion Universal Ring <b>OPTION 3:</b> ReelMatrix	<b>Pros:</b> Retention capability. <b>Cons:</b> May not retain in all circumstances. <b>Pros:</b> Similar retention capability to the Strata-G. <b>Cons:</b> May not retain in all circumstances. <b>Pros:</b> ReelMatrix will work where a sectional matrix ring will not grip appropriately. <b>Cons:</b> Less ideal contour and contact to a sectional matrix solution.
<b>Mal-Positioned</b>	<b>BEST CHOICE:</b> Quad Ring	<b>Pros:</b> Adapts to tooth rotation. <b>Cons:</b> N/A






# Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials
<b>Back-to-Back</b> 	<b>BEST CHOICE:</b> Quad Ring, Quad Wedge, Firm Band <b>Pros:</b> Ring driver tip and split wedge work together. Firm Bands prevent deformation. Reduced patient chair time. <b>Cons:</b> Learning curve.
<b>Wide Prep</b> 	<b>BEST CHOICE:</b> Strata-G Wide Prep Ring, Firm Band <b>Pros:</b> Smaller ring diameter. The only sectional ring that can perform on wide preps. <b>Cons:</b> Limited application. Not recommended for basic Class II. <b>OPTION 2:</b> 3D Fusion Wide Prep Ring, Firm Band <b>Pros:</b> Similar performance to Strata-G. <b>Cons:</b> Larger ring diameter. <b>OPTION 3:</b> Reel Matrix <b>Pros:</b> For teeth too broken down for a sectional matrix system. <b>Cons:</b> Less ideal contour and contact compared to a sectional matrix solution.
<b>Deep Prep</b> 	<b>BEST CHOICE:</b> Quad Ring, Quad Wedge, and Firm Band <b>Pros:</b> Split tip wedge sits deep in sulcus and is less likely to be pushed up by the gingiva. <b>Cons:</b> N/A <b>OPTION 2:</b> Strata-G Ring, Strata-G Wedge, and Firm Band <b>Pros:</b> Strata-G wedge sits deeper than most wedges. <b>Cons:</b> Strata-G wedge does not sit as deep as Quad.
<b>Class III</b> 	<b>BEST CHOICE:</b> Anterior Fusion <b>Pros:</b> Natural emergence profile and interproximal contours. <b>Cons:</b> Learning curve. Shape may not be appropriate for all teeth. (i.e. bulbous). Not cure through. <b>OPTION 2:</b> Vari-Strip <b>Pros:</b> Simple. Cure through. <b>Cons:</b> Susceptible to deformation during use.
<b>Class IV</b> 	<b>BEST CHOICE:</b> Anterior Fusion <b>Pros:</b> Natural emergence profile and interproximal contours. <b>Cons:</b> Learning curve. Shape may not be appropriate for all teeth. (i.e. bulbous). Not cure through. <b>OPTION 2:</b> Vari-Strip <b>Pros:</b> Simple. Cure through. <b>Cons:</b> Susceptible to deformation during use.
<b>Class V</b> 	<b>BEST CHOICE:</b> Blue View Cervical <b>Pros:</b> Compresses composite into the restoration. Limited finishing required. <b>Cons:</b> N/A

# Matrix Selection Guide

## Special Cases

Procedure	Recommended Materials
<b>Deep Margin Elevation</b> 	<p><b>BEST CHOICE:</b> ReelMatrix, Deep Margin Elevation Band (P/N RM500H)</p> <p><b>Pros:</b> Ability to seal very deep lesions. Lightweight, retainerless.</p> <p><b>Cons:</b> Will not restore whole tooth. Step one of a two step process.</p> <p><b>OPTION 2:</b> Deep Margin Elevation Tofflemire Band (P/N RM500H) with a standard tofflemire retainer</p> <p><b>Pros:</b> Ability to seal very deep lesions.</p> <p><b>Cons:</b> Retainer can torque the band out of position.</p>
<b>Last Tooth in Arch/ Missing Adjacent Tooth</b> 	<p><b>BEST CHOICE:</b> ReelMatrix</p> <p><b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. Retainerless bands allow for simultaneous multi-tooth restoration.</p> <p><b>Cons:</b> N/A</p> <p><b>OPTION 2:</b> Tofflemire Slickbands</p> <p><b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition.</p> <p><b>Cons:</b> Retainer can torque the band out of position. Cannot perform simultaneous multi-tooth restorations.</p>
<b>Core Buildup</b> 	<p><b>BEST CHOICE:</b> ReelMatrix</p> <p><b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition. Retainerless bands allow for simultaneous multi-tooth restoration.</p> <p><b>Cons:</b> N/A</p> <p><b>OPTION 2:</b> Tofflemire Slickbands</p> <p><b>Pros:</b> Circumferential band can be placed without the presence of adjacent dentition.</p> <p><b>Cons:</b> Retainer can torque the band out of position. Cannot perform simultaneous multi-tooth restorations.</p>
<b>Pedo</b> 	<p><b>BEST CHOICE:</b> Strata-G Universal Ring</p> <p><b>Pros:</b> Retention capability on shorter teeth.</p> <p><b>Cons:</b> N/A</p> <p><b>OPTION 2:</b> ReelMatrix</p> <p><b>Pros:</b> Speed.</p> <p><b>Cons:</b> Less ideal contour and contact to a sectional matrix solution.</p>
<b>Distal of Canine</b> 	<p><b>BEST CHOICE:</b> Strata-G Universal Ring</p> <p><b>Pros:</b> Retention capability.</p> <p><b>Cons:</b> May not retain in all circumstances.</p> <p><b>OPTION 2:</b> 3D Fusion Universal Ring</p> <p><b>Pros:</b> Similar retention capability to the Strata-G.</p> <p><b>Cons:</b> May not retain in all circumstances.</p> <p><b>OPTION 3:</b> ReelMatrix</p> <p><b>Pros:</b> ReelMatrix will work where a sectional matrix ring will not grip appropriately.</p> <p><b>Cons:</b> Less ideal contour and contact to a sectional matrix solution.</p>
<b>Mal-Positioned</b>	<p><b>BEST CHOICE:</b> Quad Ring</p> <p><b>Pros:</b> Adapts to tooth rotation.</p> <p><b>Cons:</b> N/A</p>