# Composi-Tight® 3D Fusion™

## Sectional Matrix System

**User instructions** (*ENGLISH*)

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**Brugsanvisninger** (*NORWEGIAN / NORGE*)

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## User instructions.

The Composi-Tight® 3D Fusion™ Sectional Matrix System is designed to create tight, anatomical contacts on Class II composite restorations. The contoured sectional matrix bands create proper tooth anatomy while three styles of rings create the required tooth separation and band adaptation for a tight, natural contact.

For Patent information see [www.garrisondental.com/patents](http://www.garrisondental.com/patents)  
Non-sterile, RX only, always use a rubber dam.

## Important Cleaning and Sterilization Instructions:

**Composi-Tight® 3D Fusion™ Soft Face™ Rings:** Remove all visible debris after use and prior to sterilization. DO NOT use solutions containing acetone or ketone. Steam autoclave sterilization only.

**Composi-Tight® 3D Fusion™ Sectional Matrix Bands:** Matrix bands are single use disposable items. Sterilization prior to use is recommended and can be performed using steam autoclave.

**Instruments & Forceps:** Instruments may be dry heat sterilized, steam autoclaved or cold sterilized.

**Composi-Tight® 3D Fusion™ Wedges:** Wedges are single use, disposable items. Spray disinfect only. Do not autoclave.

**For complete sterilization instructions go to: [garrisondental.com/sterilization](http://garrisondental.com/sterilization)**

## Technique for Class II restorations utilizing Composi-Tight® 3D Fusion™ Soft Face™ Rings

- Place rubber dam.** Isolation of the area helps to prevent contamination from blood and saliva. As with many dental procedures, aspiration of small parts is a valid concern.
- Place FenderWedge®** prior to preparation to pre-separate and protect the adjacent tooth. (figure 1) Prepare tooth with conservative Class II preparation. Once the preparation is complete, remove the FenderWedge® and place your matrix.
- Place Composi-Tight® 3D Fusion™ sectional matrix band.** Select a matrix band that most closely approximates the occlusogingival height of the tooth. The Grab-Tab should be oriented toward the occlusal margin and may be bent over the adjacent tooth after placement. (figure 2)
- Insert Composi-Tight® 3D Fusion™ wedge.** Keep a finger on the matrix band to prevent the wedge from dislodging the band during insertion. Wedge firmly, a high viscosity resin packed into the preparation may displace the band away from the tooth. (figure 3)
- Apply the Composi-Tight® 3D Fusion™ Soft Face™ ring retainer.** Select a retainer ring that is appropriate for the embrasure and tooth being restored. Spread it with the \*Composi-Tight 3D Fusion™ Ring Placement forceps and place it over the wedge. (figure 4)
- Very Important:** push the matrix band firmly into contact with the adjacent tooth in the desired contact area.
- Use your preferred etching, disinfecting and bonding technique** following manufacturer's instructions. Apply your preferred composite resin following manufacturer's instructions.
- Remove the ring, wedge and band.**
- Contour and polish the restoration.**


figure 1

figure 2


figure 3


figure 4


figure 5

## Stacking of Rings & MOD Restorations

The orange tipped and blue tipped versions of the Composi-Tight® 3D Fusion™ Soft-Face™ Rings are angled to allow stacking of rings for performing M.O.D.s and multiple-tooth restorations. Rings may be stacked over each other. (figure 5)

\*(Garrison part number FXP01 Composi-Tight® 3D Fusion™ Ring Placement Forceps will hold the rings more securely and provide ample force to open the rings. NOTE: Most rubber dam clamp forceps will not open the Composi-Tight® 3D Fusion™ Rings far enough to allow placement.)

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