

Cervical Decay, Ionomer Filling and Matrix System: An Alternative to Composite

Procedure/Study by

Dr. David Gerdolle, Switzerland

Composi-Tight® 3D XR 



Cervical Decay, Ionomer Filling and Matrix System: An Alternative to Composite

By Dr. David Gerdolle, Switzerland

David Gerdolle was born in Epinal (France) in May 1970, and graduated from the Dental Faculty of Nancy in 1993. Since then, he has achieved a number of PGs in Oral and Osteo-articular Biology, Prosthodontics and Adhesive Dentistry.

David was a member of the teaching staff of the University of Nancy from 1993 until 2005, and has been in private practice in Vevey-Montreux (CH) since 2006. His practice is dedicated to conservative and minimally invasive dentistry.

David is involved in postgraduate dental teaching at the University of Paris (F), as well as in research programs on resin composites at the University of Nancy. He is lecturing in seminars and congresses, giving hands-on courses and performing live case demonstrations in the field of adhesive dentistry and indirect restorations in particular.


As author or co-author, David has published widely in scientific national and international journals on this topic.

01



Pre-op X-ray: cervical decay located below the enamel-cementum junction

02



Cavity ready to be filled: Garrison Composi-Tight® matrix system in place

03



Due to the location of the decay and to the poor interproximal hygiene of the patient, it has been decided to fill the cavity using a glass ionomer cement (Equia Forte from GC). It has been previously checked that the distal part of the treated tooth was not functional.

04



After three minutes, the glass ionomer has set. The matrix is then removed and the material is finished and coated (Equia Forte Coat, GC).

05



Post-op view

06



Post-op X-ray



Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456
Toll free 888.437.0032

European office:
Carlstr. 50
52531 Uebach-Palenberg
Germany
Tel. +49 (0) 2451 971 409

www.garrisondental.com

CS-000-02-93