



Garrison
Dental Solutions

2021 PRODUCT CATALOG

table of contents

MATRIX SYSTEMS

sectional matrix systems	Composi-Tight® 3D Fusion™	4-5
	NiTin™	6-7
	Other sectional systems and refills	8-9
retainerless Tofflemire-style	Reel Matrix™	30-31
	Pinch™ Matrix	34
traditional Tofflemire-style	Slick Bands™	32-33
anterior matrix systems	VariStrip™	36
cervical matrix systems	Blue View™ Cervical Matrices.	37

WEDGES

plastic	Composi-Tight® 3D Fusion™ Wedge	14
	A+Wedge®	15
	WedgeWands®	15
	G-Wedge™	16
wood	SoftWedge™	16
pre-wedge & shield	FenderWedge®	17

IMPORTANT ACCESSORIES

polishing systems	Rally™	24-25
	eZr™	26-27
interproximal finishers	FitStrip™	20-23
forceps	Composi-Tight® 3D Fusion™ Forceps	39
	Composi-Tight® Matrix Forceps	38
	NiTin™ Forceps	38
	Ring Placement Forceps	39
	Universal Placement Forceps	39
specialty	PerForm™ Contact Former	35
composite instruments	Multi-Function Non-Stick	40
	Standard Non-Stick	41


INFORMATION PAGES

Composi-Tight® 3D Fusion™	2-3
Band Comparison Chart	10-12
Matrix Selection Guide	42-43
Order Form	44
Ordering Information	45
Guarantee	back cover

Only a sectional matrix can deliver results like this!

Tofflemire-style System

- ❌ Fails to restore proximal anatomy
- ❌ Thin contact at the marginal ridge
- ❌ Large food trap below
- ❌ Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Composit-Tight® 3D Fusion™ Sectional Matrix System

- ✓ Operator-friendly retaining system
- ✓ Naturally contoured bands
- ✓ Anatomically correct contacts
- ✓ Contacts at the height of contour
- ✓ Rings produce optimal tooth separation for excellent, tight contacts

Only Composit-Tight® 3D Fusion™ can make the technique this intuitive and frustration free!



Place the matrix band:

Dead-soft three dimensionally contoured matrices have placement tabs for easier occlusal placement and Garrison's revolutionary non-stick coating for hassle-free removal.



Wedge firmly:

Composit-Tight® 3D Fusion™ Ultra-Adaptive wedges glide in but won't back out. Their soft, flexible exterior readily adapts to interproximal irregularities for a superior gingival seal.



Place the ring:

All three Composit-Tight® 3D Fusion™ separator rings will fit right over the wedge for easier placement. With Ultra-Grip™ retention extensions, they will stay firmly in place from the distal of the canine to the widest of molar preparations.

Composi-Tight® 3D Fusion™ Cases

Great adaptation to the teeth eliminates the need for adjustment with polishers in the proximal areas.



Images courtesy of Clarence Tam, Auckland, NZ

Wide preps can be handled in one simple step when you have the right (FX600 green) ring.



Images courtesy of Dr. Nikhil Sethi, BDS (Lon) Msc (Lon), Birmingham, UK

Back-to-back requires a shift in technique—we can show you how.*



Images courtesy of Dr. Jeffrey Kelley, DDS, Zeeland, MI

Quadrants can be a challenge. Ease your concerns with tips and tricks from the pros.



Images courtesy of Dr. Mike Malone, DDS, Lafayette, LA

*To see full writeups, cases and reviews go to www.garrisondental.com/learning-center
Do you have a case you want to share? Email it to us at contact@garrisondental.com.

The Award Winning Composi-Tight® 3D Fusion™ Sectional Matrix System



TESTIMONIAL

"3D rings are the real topic of Garrison's systems. The "v" shape of the ring that fits in the interproximal area allows a good fitting between the cavity margins and the matrix in the buccal and palatal walls. This results in an easier positioning of the composite masses close to the cavity margins and also the final remodeling usually necessary at the time of removal of the matrix will be very minimal. The rings also permit a divergence of the interproximal dental elements, and this causes a great point of contact. Garrison systems make Class II restorations simpler and more predictable and also reduce the operating time of the finishes when the matrix is taken off."

Enrico Cogo, DDS
Legnago, Ferrara, and Goito, Italy

Composi-Tight® 3D Fusion™

**A system designed
to handle even the
toughest cases**

More than 20 years of sectional matrix leadership and an incredible amount of feedback from users all over the globe were poured into the design. All this hard work has resulted in separator rings with significantly improved performance and ease of use.



Hugs the curves – Soft-Face™ silicone tips with enhanced marginal ridge anatomy hug the matrix band to the tooth for flash free restorations

It just keeps going... Enhanced durability for extended ring life and greater tooth separation

Wider indications, wider smiles!



Dentistry's first Wide Preparation separator ring (green) makes Composi-Tight® 3D Fusion™ the most versatile and user-friendly sectional matrix system available. The Wide Prep ring greatly simplifies what until now has been one of the most challenging posterior composite restorations. Now you can restore these big preps with confidence knowing that you will achieve excellent contact and contour without having to take extra steps.

Matrix band perfection

Composi-Tight® 3D Fusion™ Full Curve Non-Stick Matrices are three dimensionally contoured, ultra-thin and dead-soft to help you recreate ideal interproximal anatomy.



Easy in – easy out! A Grab-Tab™ centered on the occlusal edge of the band simplifies placement while Slick Bands™ non-stick coating reduces composite adhesion by an impressive 92%!

Go deep or

wide – Improved subgingival extensions and added wrap-around length make for the perfect matrix band to use with the green Wide Preparation ring. Now you're ready to tackle even the widest preps!



FXR-KFF-11

See a product demo online!



Composi-Tight® 3D Fusion™

Sectional Matrix System



Composi-Tight® 3D Fusion™ Bands # .0016" thick	FX100	FX150	FX175	FX200	FX300
Non-stick dead soft stainless steel	bicuspid 4.0 mm w/ext 4.4 mm	bicuspid with deep cervical preparations 4.1 mm w/ext 6.0 mm	bicuspid and small molars 5.0 mm w/ext 5.6 mm	molars 5.4 mm w/ext 6.6 mm	molars with deep cervical preparations 6.0 mm w/ext 8.7 mm

Kit

- 1 Short Ring (blue)
- 1 Tall Ring (orange)
- 1 Wide Prep Ring (green)
- 70 Assorted Full Curve Non-Stick Matrix Bands
- 80 Assorted Ultra-Adaptive Wedges
- 4 Assorted FenderWedges® (pg 17)
- 12 Assorted Rally™ Polishers (pg 24–25)
- 1 Ring Placement Forceps (pg 39)
- 1 BFA-1 Band Placement Forceps (pg 38)
- 1 TN009 Multi-Function Instrument (pg 40)

Item

Item	Price
FXR-KFF-11	\$779
FXR-KFF-10 (BFA-1 and TN009 not in kit)	\$569
FXR-KFF-14 (TN009 not in kit)	\$689

Composi-Tight® 3D Fusion™ Rings

- Short Ring (blue)
- Tall Ring (orange)
- Wide Prep Ring (green)
- 3-Ring Kit (1 of each ring)

Item	Qty.	Price
FX400	2	\$235
FX500	2	\$235
FX600	2	\$235
FXR01	3	\$344

Composi-Tight® 3D Fusion™ Bands

- Bicuspid
- Bicuspid w/ deep cervical preps
- Bicuspid and small molars
- Molars
- Molars w/ deep cervical preps

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
Bicuspid	50	\$81	100	\$142
Bicuspid w/ deep cervical preps	30	\$81	60	\$142
Bicuspid and small molars	50	\$81	100	\$142
Molars	50	\$81	100	\$142
Molars w/ deep cervical preps	30	\$81	60	\$142

NEW!

NiTiⁱⁿ

NiTinol Sectional Matrix System

**NiTinol + PEEK
for the ultimate
nickel titanium
ring performance**

As the name implies, the new NiTin[™] sectional matrix system separator rings are made from a special formulation of NiTinol nickel titanium for increased strength and resiliency. Rather than typical material used by other brands, NiTin[™] is made with drawn-wire Nitinol which produces superior spring characteristics. Add in the PEEK reinforcement on the bow of the ring and you've got a system that produces exceptional contacts for hundreds and hundreds of uses.


Outperforms all other higher priced competitors' rings!



Smaller, less "bulky" design fits well in tight spots and is still stackable for MODs.

Drawn-wire nickel titanium for the ultimate in long life!


Proprietary heat-set processing increases strength and resiliency.



PEEK (polyether ether ketone) super-polymer reinforcement reduces ring stretch.

Ultra-durable PEEK tips are the strongest, longest lasting polymer tips on any ring.

Precisely engineered tip geometry allows NiTin[™] to grip where others pop-off—all while reducing flash.





NiTinol Sectional Matrix System



NiTIn™ Bands # .0016" thick	NTG100	NTG150	NTG175	NTG200	NTG300
Non-stick dead soft stainless steel	bicuspid 4.0 mm w/ext 4.4 mm	bicuspid with deep cervical preparations 4.1 mm w/ext 6.0 mm	bicuspid and small molars 5.0 mm w/ext 5.6 mm	molars 5.4 mm w/ext 6.6 mm	molars with deep cervical preparations 6.0 mm w/ext 8.7 mm

Kit

- 1 Long NiTin™ Ring
- 1 Standard NiTin™ Ring
- 100 Assorted NiTin™ Matrices
- 100 Assorted G-Wedges™
- 1 NiTin™ Ring Placement Forceps
- 1 NiTin™ Band Forceps

Item

NTG-KMG-01

Price

\$455

NiTIn™ Rings

- Long Ring (black)
- Standard Ring (white)

Item

NTG400
NTG500

Qty.

2
2

Price

\$205
\$205

NiTIn™ Bands

- Bicuspid
- Bicuspid w/ deep cervical preps
- Bicuspid and small molars
- Molars
- Molars w/ deep cervical preps

Item

NTG100
NTG150
NTG175
NTG200
NTG300

Qty.

50
50
50
50
50

Price

\$64
\$64
\$64
\$64
\$64

The Complete NiTin™ System

NTG-KMG-01



Composi-Tight® 3D XR

Sectional Matrix System



Slick Bands™ XR .0016" thick	SXR100	SXR150	SXR175	SXR200	SXR300
Non-stick dead soft stainless steel	bicuspids 4.6 mm	bicuspids with deep cervical preparations 3.8 mm w/ext 5.4 mm	bicuspids and small molars 5.4 mm	molars 6.4 mm	molars with deep cervical preparations 6.4 mm w/ext 8.6 mm

Kit

3D XR System Kit

- 1 Composi-Tight® 3D XR Ring
- 2 Original 3D500 Soft-Face™ Rings
- 100 Slick Bands™ XR Non-stick Matrices
(30 ea. SXR100, SXR200, SXR175; 5 ea. SXR150, SXR300)
- 100 WedgeWands® (25 ea. size)
- 20 FenderWedges® (5 ea. size)
- 1 Ring Placement Forceps (AUMRDF-100))

Item

Price

3D-KSW-10 \$569

3D Rings

Soft-Face™ Extra Retention Rings
Soft-Face™ Rings

Item

Qty.

Price

3DXR 2 \$229
3D500 2 \$224

Composi-Tight® SXR Bands	Item	Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Bicuspids	SXR100	50	\$73	100	\$127
Bicuspids w/ deep cervical preps	SXR150	30	\$73	60	\$127
Bicuspids and small molars	SXR175	50	\$73	100	\$127
Molars	SXR200	50	\$73	100	\$127
Molars w/ deep cervical preps	SXR300	30	\$73	60	\$127

Composi-Tight® 3D Clear

Cure-through sectional matrix system



Composi-Tight® 3D Clear .0020" thick	MC100	MC175	MC200
Blue tinted polyester	premolar 5.0 mm	bicuspids and small molars 5.9 mm	tall molar 6.8 mm

Kit

3D CLEAR System Kit

- 1 Composi-Tight® 3D XR Ring
- 2 3D600 Clear Soft-Face™ Rings
- 100 Cure-Through matrices
(30 ea. MC100, MC200; 40 ea. MC175)
- 100 Cure-Through WedgeWands® (25 ea. size)
- 20 FenderWedges® (5 ea. size)
- 1 Ring Placement Forceps (AUMRDF-100)

Item

Price

3D-KCC-10 \$569

3D Clear Rings

Soft-Face™ Extra Retention Rings
Clear Soft-Face™ Rings

Item

Qty.

Price

3DXR 2 \$229
3D600 2 \$224

3D Clear Bands	Item	Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Premolar matrix bands	MC100	50	\$73	100	\$127
Bicuspid and small molar matrix bands	MC175	50	\$73	100	\$127
Tall molar matrix bands	MC200	50	\$73	100	\$127

To order standard packs, add a -M after the item #.
Complete band selection guide pages 10, 11 and 12.

Composi-TightGold®

Sectional Matrix System



Composi-Tight Gold # .0015" thick	AU050	AU100	AU150	AU200	AU300
Gold colored dead soft stainless steel	pediatric 3.2 mm w/ext 4.2 mm	bicuspid 4.6 mm	bicuspid with deep cervical preps 3.8 mm w/ext 5.4 mm	molars 6.4 mm	molars with deep cervical preps 6.4 mm w/ext 8.6 mm

Kit

Gold System w/Forceps

- 250 Matrix Bands (100 ea. AU100, AU200; 25 ea. AU150, AU300)
- 6 G-Rings (3 standard, 3 long)
- 1 Ring Placement Forceps

Item

AUK2

Price

\$490

Gold G-Rings™

G-Rings™ w/Standard Length Tines

G-Rings™ w/Long Length Tines

Item

AU400

Qty.

3

Price

\$98

AU500

3

\$98

Gold Bands

	Item	Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Pediatric	AU050	25	\$60	50	\$92
Bicuspid	AU100	50	\$60	100	\$92
Bicuspid w/ deep cervical preps	AU150	25	\$60	50	\$92
Molars	AU200	50	\$60	100	\$92
Molars w/ deep cervical preps	AU300	25	\$60	50	\$92

Forceps

Ring Placement Forceps



Item

AUMRDF-100

Price

\$83

Composi-Tight™

Sectional Matrix System



Composi-Tight # .0015" thick	B050	B100	B150	B200	B300
Polished dead soft stainless steel	pediatric 3.2 mm w/ext 4.2 mm	bicuspid 4.6 mm	bicuspid with deep cervical preps 3.8 mm w/ext 5.4 mm	molars 6.4 mm	molars with deep cervical preps 6.4 mm w/ext 8.6 mm

Kit

Original System w/Forceps

- 225 Matrix Bands (100 ea. B100, B200; 25 B300)
- 7 G-Rings (4 standard, 3 long)
- 1 Ring Placement Forceps

Item

K102

Price

\$450

Original G-Rings™

G-Rings™ w/Standard Length Tines

G-Rings™ w/Long Length Tines

Item

G100

Qty.

3

Price

\$88

G200

3

\$88

Original Bands

	Item	Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Pediatric	B050	25	\$52	50	\$79
Bicuspid	B100	50	\$52	100	\$79
Bicuspid w/ deep cervical preps	B150	25	\$52	50	\$79
Molars	B200	50	\$52	100	\$79
Molars w/ deep cervical preps	B300	25	\$52	50	\$79

To order standard packs, add a -M after the item #.

Complete band selection guide pages 10, 11 and 12.

Composi-Tight® 3D Fusion™

NEW!

Firm Matrix Bands

New **Composi-Tight® 3D Fusion™ Firm** matrix bands with **eZ-Place™** technology maintain their shape during placement making them *significantly easier* to insert!

eZ-Place™ technology combines specially processed stainless steel with a proprietary matrix band forming process to yield Fusion Firm's remarkable handling characteristics. They resist deformation nearly 70% better than typical dead-soft bands yet are pre-contoured and super thin at 0.0015". You're ensured the tight contacts and ideal anatomy you've come to expect from the Composi-Tight® 3D Fusion™ system.

Fusion Firm:

- Nearly 70% better at resisting deformation!
- Super-thin at only 0.0015"
- Five ideal sizes!

Perfect for:

- Tight interproximal spaces!
- More conservative preps!
- First-time users and veterans alike!

Description	Item	Qty.	Kit
Assorted Band Kit: 5 sizes	FXHB04	300	\$355
Assorted Band Mini Kit: 5 sizes	FXHB05	150	\$185



FXH Series

Composi-Tight® 3D Fusion™ Firm						
	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
thickness 0.0015"						
height						
NA						
		FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
FXH100	50	\$78	100	\$137
FXH150	50	\$78	100	\$137
FXH175	50	\$78	100	\$137
FXH200	50	\$78	100	\$137
FXH300	50	\$78	100	\$137

FX Series

Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.						
	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
thickness 0.0016"						
height						
NA						
		FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm


Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
FX100	50	\$81	100	\$142
FX150	30	\$81	60	\$142
FX175	50	\$81	100	\$142
FX200	50	\$81	100	\$142
FX300	30	\$81	60	\$142

FX Series

SXR Series

B Series

M – SM – AU Series



Band Comparison Chart

NTG Series

Pediatric | Bicuspid | Bicuspid with deep cervical preps | Bicuspid and small molars | Molars | Molars with deep cervical preps

NiTi™ NiTi™ Full contour matrix bands for ideal anatomy.

thickness
0.0016"

height

NA



NTG100
4.0 mm
w/ext 4.4 mm



NTG150
4.1 mm
w/ext 6.0 mm



NTG175
5.0 mm
w/ext 5.6 mm



NTG200
5.4 mm
w/ext 6.6 mm



NTG300
6.0 mm
w/ext 8.7 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
NTG100	25	\$36	50	\$64
NTG150	25	\$36	50	\$64
NTG175	25	\$36	50	\$64
NTG200	25	\$36	50	\$64
NTG300	25	\$36	50	\$64

SXR Series

Slick Bands™ XR Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.

thickness
0.0016"

height

NA



SXR100
4.6 mm



SXR150
3.8 mm
w/ext 5.4 mm



SXR175
5.4 mm



SXR200
6.4 mm



SXR300
6.4 mm
w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
SXR100	50	\$73	100	\$127
SXR150	30	\$73	60	\$127
SXR175	50	\$73	100	\$127
SXR200	50	\$73	100	\$127
SXR300	30	\$73	60	\$127

SM Series

Slick Bands™ Slick Bands™ Matrices have non-stick coating.

thickness
0.0016"

height



SM050
3.2 mm
w/ext 4.2 mm



SM100
4.6 mm



SM150
3.8 mm
w/ext 5.4 mm



SM175
5.5 mm



SM200
6.4 mm



SM300
6.4 mm
w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
SM050			60	\$127
SM100	50	\$73	100	\$127
SM150	30	\$73	60	\$127
SM175	50	\$73	100	\$127
SM200	50	\$73	100	\$127
SM300	30	\$73	60	\$127







To order standard packs, add a -M after the item #.

AU Series

	Pediatric	Bicuspsids	Bicuspsids with deep cervical preps	Bicuspsids and small molars	Molars	Molars with deep cervical preps
Composi-Tight Gold®	The Gold bands are the same as the Composi-Tight® bands.					
thickness 0.0015"				NA		
height	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm		AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
AU050	25	\$60	50	\$92
AU100	50	\$60	100	\$92
AU150	25	\$60	50	\$92
AU200	50	\$60	100	\$92
AU300	25	\$60	50	\$92

M Series

Composi-Tight®	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.					
thickness 0.0015"						
height	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm




Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
M050	25	\$60	50	\$92
M100	50	\$60	100	\$92
M150	25	\$60	50	\$92
M175	50	\$60	100	\$92
M200	50	\$60	100	\$92
M300	25	\$60	50	\$92

B Series

Composi-Tight® Original	The original bands have a slightly shorter length.					
thickness 0.0015"				NA		
height	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm		B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
B050	25	\$52	50	\$79
B100	50	\$52	100	\$79
B150	25	\$52	50	\$79
B200	50	\$52	100	\$79
B300	25	\$52	50	\$79

MC Series

Composi-Tight® 3D Clear	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.					
thickness 0.0020"						
height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA

Item	Standard (-M)		Bulk	
	Qty.	Price	Qty.	Price
MC100	50	\$73	100	\$127
MC175	50	\$73	100	\$127
MC200	50	\$73	100	\$127

To order standard packs,
add a -M after the item #.

Cases using Garrison® Wedges



Images courtesy of Dr. Mike Malone, DDS, Lafayette, LA



Images courtesy of Dr. Robert Rosenfeld, DDS, Westwood, NJ

Composit-Tight® 3D Fusion™ Wedge

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.



Images courtesy of Dr. Robert Rosenfeld, DDS, Westwood, NJ

Wedge Wands®

You can twist the handle off or leave it on. If you choose to leave the handle on, it can be used to snug the wedge in tighter once the ring separates the teeth. Then remove it when finished without ever having to pass an instrument back and forth with your assistant.



Images courtesy of Dr. Schuyler Van Gorden, DDS, Eau Clair, WI



Images courtesy of Dr. Strassler, DDS

FenderWedge®

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations. FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and pre-separates the teeth, thus simplifying the application of a matrix.



Images courtesy of Dr. David S. Hornbrook, DDS, FAACD, San Diego, CA

TESTIMONIALS

"Garrison's redesigned **3D Fusion™ Wedges** feature soft "fins" that smoothly adapt matrix bands to tooth irregularities and prevent the wedge from backing out. These are helpful in a wide array of situations, to eliminate the likelihood of flash at the gingival margin."

*Dr. Robert Rosenfeld, DDS
Westwood, NJ*

"The Garrison **A+Wedge™** minimized or eliminated bleeding. It's particularly beneficial on patients with slight to mild gingivitis."

*Joe Arbutante, DDS
Garrison Advisory Board Member*


"Even those with steady hands and stellar technique occasionally and inadvertently nick the proximal tooth or restoration as they are preparing the adjacent tooth. I know I have! The **FenderWedge®** is a simple device that I wish I would have had the opportunity to use since the day I started practicing dentistry. It would have saved me much aggravation and time. It is now a "must-have" item during all my restorative procedures."

*David S. Hornbrook, DDS, FAACD
International Lecturer and
Private Practice, San Diego, CA*

Composi-Tight® 3D Fusion™

Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



A real softie Composi-Tight® 3D Fusion's Soft-Face™ over-mold allows the wedge to do what no other wedge can truly do – actually adapt to interproximal irregularities.

Tough on the inside
The firm inner skeleton allows for easy insertion and tooth separation just like a traditional wedge.



Fins with serious function
Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.



Composi-Tight® 3D Fusion™ Wedge

Description

Composi-Tight® 3D Fusion™ Wedge Kit
Composi-Tight® 3D Fusion™ Wedge Kit

Item	Qty.	Kit
FXK4	400	\$207
FXK4-M	200	\$109

Refills

Item	Standard		Item	Bulk		
	Qty.	Price		Qty.	Price	
X-Small Yellow	FXYL-M	50	\$32	FXYL	100	\$59
Small Blue	FXBL-M	50	\$32	FXBL	100	\$59
Medium Orange	FXOR-M	50	\$32	FXOR	100	\$59
Large Green	FXGR-M	50	\$32	FXGR	100	\$59



- Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva
- Improves clarity of restoration area
- Eliminates need to apply an astringent separately

A+Wedge®

Description	Item	Qty.	Kit
A+Wedge® Kit	GWAK4	200	\$114

Refills	Item	Bulk	
		Qty.	Price
X-Small Yellow	GWAYL	100	\$59
Small Blue	GWABL	100	\$59
Medium Orange	GWAOR	100	\$59
Large Green	GWAGR	100	\$59

Insert



Wait & Restore



Remove & Rinse



WedgeWands®



WedgeWands
also available in
Cure-Through

WedgeWands®

Description	Opaque / Clear	Qty.	Kit
WedgeWands® Kit	WK4 / WCK4	400	\$178

Refills	Opaque / Clear	Standard		Opaque / Clear	Bulk	
		Qty.	Price		Qty.	Price
X-Small Yellow	SWYL / SWCYL	100	\$53	WYL / WCYL	300	\$143
Small Blue	SWBL / SWCBL	100	\$53	WBL / WCBL	300	\$143
Medium Orange	SWOR / SWCOR	100	\$53	WOR / WCOR	300	\$143
Large Green	SWGR / SWCGR	100	\$53	WGR / WCGR	300	\$143



G-Wedge™

Perfect wedge. Perfect contact.

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.

G-Wedge™

Description	Item	Qty.	Kit
G-Wedge™ Kit	GWK4	400	\$153

Refills	Item	Standard		Item	Bulk	
		Qty.	Price		Qty.	Price
X-Small Yellow	GSWYL	100	\$45	GWYL	300	\$121
Small Blue	GSWBL	100	\$45	GWBL	300	\$121
Medium Orange	GSWOR	100	\$45	GWOR	300	\$121
Large Green	GSWGR	100	\$45	GWGR	300	\$121

SoftWedge™

Superior adaptation and seal



The first wooden wedge designed for use with a sectional matrix system.

Soft Wedges are wooden interproximal wedges that have been optimized for use with modern restorative materials. The ability to compress during insertion, then swell to help fill voids ensures that your preparations stay clean and dry—an absolute must for today's bonded restorations!

Soft Wedge™

Description	Item	Qty.	Kit
Soft Wedge™ Kit	WDK4	400	\$95

Refills	Item	Standard		Item	Bulk	
		Qty.	Price		Qty.	Price
X-Small Yellow	WSDYL	100	\$29	WDYL	300	\$79
Small Blue	WSDBL	100	\$29	WDBL	300	\$79
Medium Orange	WSDOR	100	\$29	WDOR	300	\$79
Large Green	WSDGR	100	\$29	WDGR	300	\$79

FENDERWEDGE®

Interproximal Tooth Shields

Pre-separate and Protect

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations.

FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and separates the teeth, thus simplifying the application of a matrix.



"Even those with steady hands and stellar technique occasionally and inadvertently nick the proximal tooth or restoration as they are preparing the adjacent tooth. I know I have! The Fender Wedge is a simple device that I wish I would have had the opportunity to use since the day I started practicing dentistry. It would have saved me much aggravation and time. It is now a "must-have" item during all my restorative procedures."

David S. Hornbrook, DDS, FAACD
International Lecturer and
Private Practice, San Diego, CA

FenderWedge®

Description	Item	Qty.	Kit
FenderWedge® Kit	TSK01	80	\$128

Refills	Item	Standard		Item	Bulk	
		Qty.	Price		Qty.	Price
X-Small Purple	TSPR-M	30	\$56	TSPR	60	\$102
Small Orange	TSOR-M	30	\$56	TSOR	60	\$102
Medium Green	TSGR-M	30	\$56	TSGR	60	\$102
Large Yellow	TSYL-M	30	\$56	TSYL	60	\$102



FenderWedge is NOT a sectional matrix system nor does it replace the use of a matrix band.

Great technique tips and product help are only a click away!

Join thousands of subscribers on our YouTube channel and get notified instantly when new content is posted.

The Garrison® YouTube channel is NOT a collection of infomercials. It's all about helpful information on a wide variety of topics. From short two-minute tips on how to reduce flash on Class IIs, to more in depth live patient case presentations, the Garrison® channel delivers great content.

Go to youtube.com/garrisondental or search Garrison Dental on youtube and subscribe today!



Speaking of great content...

Get CE credits direct from Garrison®!

Go to the Learning Center at garrisondental.com to see a list of pre-recorded and upcoming live webinars – many with available CE.

Register your email (big red box in the bottom right-hand corner) to get notifications of new CE events delivered right to your inbox!

Finishing and Polishing Cases

FitStrip™ case



Images courtesy of Dr. Jeffery Rosenberg, Philadelphia, PA

eZr™ case



Rally™ case



TESTIMONIAL

"Interproximal finishing and polishing is a 'must' whenever restoring a Class II, III or IV. However, it is often difficult to achieve. Garrison's new FitStrip™ device is ideal for interproximal trimming (even light orthodontic IPR) and composite finishing."

*Dr. Howard S. Glazer, DDS FAGD
Fort Lee, NJ*

"I use the eZr™ system on all zirconia and e.max crowns. This product has been fantastic at polishing zirconia...much better than many other brands!"

*Dr. Phillip Howard, DDS
Seattle, WA*



Flexible
Interproximal finishers with a
Twist


New

Subgingival Trimmers and Finishers!

Subgingival FitStrips are perfectly designed to remove composite overhangs and excess cement quickly and safely from hard-to-reach subgingival areas. These strips are 33% more narrow at only 2.5 mm in height versus standard FitStrips, and with their center, non-abrasive Safety Zone, slip easily between teeth and below the contact area. The medium, fine and super fine durable diamond coating makes quick work of the cleanup and provides a silky-smooth finish.

Do the Twist! Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place.

One handy handle! FitStrip's simple, attachable handle makes a world of difference for both clinician ease-of-use and patient comfort.



Stay in the zone

- Non-abrasive Safety Zone protects contact area
- Super-thin – only 0.05 mm!
- Slips easily and safely between teeth

Narrow focus

- 33% more narrow than standard FitStrips
- Only 2.5 mm in height

Ready for a marathon! FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

Starts bold — finishes clean

1. Medium grit for fast removal
2. Fine grit smooths things out
3. Super fine strips put the perfect finishing touch on things

FitStrip™

Flexible
Interproximal finishers with a
Twist

Simple Subgingival Technique:

1. Slide between teeth. Starting with the medium grit (blue), twist the colored barrel until the strip is locked straight. Position the Safety Zone in the affected interproximal space. With gentle apical pressure and a slight back-and-forth motion, slide the strip between the teeth*. The barrel may be twisted into a curved position at this time to allow the strip to follow the natural tooth contours.
2. Remove the overhang/excess cement. Move Fitstrip™ in a gentle buccal-lingual back-and-forth motion. Remove FitStrip™ and test the area with floss to ensure the desired result. Repeat as necessary.
3. Finish and polish. Place the fine grit FitStrip™ (red) into the interproximal space using the method in step one. Move it in a gentle buccal-lingual back-and-forth motion three or four times. Remove the fine strip and repeat the process with the super fine (yellow) for an ideal polish.

*If the teeth are bonded together or the contact is too tight to get the 0.05 mm Safety Zone between the teeth, use the FitStrip™ serrated saw (gray barrel) first to gain access.

NEW!

FitStrip™ Subgingival Finishing and Single Sided Kit

FPSK08 contains:

2 handles

10 Strips:

2 subgingival super fine

2 subgingival fine

1 subgingival medium

1 serrated saw

1 each single sided

.08 mm, .10 mm, .13 mm, .18 mm

FitStrip™

Kits All kits come in a package of 10 FitStrips with 2 handles.

NEW Subgingival Finishing Kit

3 subgingival super fine, 3 subgingival fine, 3 subgingival medium, 1 serrated saw

NEW Subgingival Finishing and Single Sided Kit

2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)



FPSK08

NEW Subgingival Finishing Kit 3 subgingival super fine, 3 subgingival fine, 3 subgingival medium, 1 serrated saw	FPSK07	\$179
NEW Subgingival Finishing and Single Sided Kit 2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	FPSK08	\$179

FitStrip™

Flexible
Interproximal finishers with a
Twist



FPSK05

IPR Kit

FPSK05 contains:

2 handles

1 IPR gauge

10 Double sided strips (2 of each):

11 mm, .21 mm, .30 mm, 40 mm, .50 mm

Your Complete Solution for Interproximal Reduction (IPR)!

.40mm and .50mm strips easily accommodate maximum reduction cases

IPR gauges for exact reduction

IMPROVED! Attachable handle for greater ease of use and patient comfort

FitStrip™

Kits All kits come in a package of 10 FitStrips with 2 handles.

Universal Kit for IPR and finishing/polishing

1 IPR gauge

4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK04 \$197

IPR Kit

1 IPR gauge

10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK05 \$197

Starter Kit for finishing/polishing

2 serrated

4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)

FPSK01 \$179

Single Sided Kit

2 serrated

8 single sided (2 of each .08 mm, .10 mm, .13 mm, .18 mm)

FPSK02 \$179

Double Sided Kit

2 serrated

8 double sided (2 of each .11 mm, .15 mm, .21 mm, .30 mm)











FPSK03 \$179



FitStrips — a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**



Refills

		Single Sided	Double Sided	10 pk w/2 handles	4 ct packs (add – 4 to model #)
Super fine 15 Micron		FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick	Either Refill \$179	\$85
Fine 30 Micron		FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick	Either Refill \$179	\$85
Medium 46 Micron		FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick	Either Refill \$179	\$85
Coarse 90 Micron		FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick	Either Refill \$179	\$85
Serrated		FPSSAW 0.05 mm thick		Refill \$179	\$85
Medium 46 Micron			FPSTMDS 0.40 mm thick	Refill \$179	\$85
Medium 46 Micron			FPSXTMDS 0.50 mm thick	Refill \$179	\$85
Super fine 15 Micron		FPSGXFSS 0.08 mm thick	N/A	Refill \$179	\$85
Fine 30 Micron		FPSGFSS 0.10 mm thick	N/A	Refill \$179	\$85
Medium 46 Micron		FPSGMSS 0.13 mm thick	N/A	Refill \$179	\$85
FitStrip replacement handles (2-pack) FPSHANDLE-2				\$30	
IPR Gauge FPSGAUGE				\$40	



Mini Polishers for Composite

Rally gets you to the finish line faster!

The composite finish line that is. One lightning fast step will put a beautiful satin finish on posterior restorations or add the gray high shine polisher for a lustrous gloss on facials and anteriors.

And if you want to be quick off the starting line, use the blue coarse grit for rapid reduction and shaping. Rally Mini Polishers get you across the finish line and on to your next patient in record time!



Get to the point – Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.



A clean get-a-way!



You don't want to leave a trail of crumbs or spattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.



Mini and mighty!


Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.






**Remove –
Choose the Blue
(silicon carbide)**

Polishers for fast
removal and
contouring



**Smooth –
Choose the Red
(aluminum oxide)**

Polishers for
smoothing and silky
shine



**Shine –
Choose the
Gray (diamond
impregnated)**

Polishers for ultra
high shine



*Slow speed right
angle latch shank.
Use with water.*

Rally™

Item	Description	Qty.	Price
FPCCK01	Composite Polisher Complete Kit	60 (10 each)	\$151


Item	Refills	Qty.	Price
FPCCK010	Blue coarse grit point composite polisher	30	\$59
FPCCK040	Blue coarse grit cup composite polisher	30	\$59
FPCCK010	Red medium grit point composite polisher	30	\$59
FPCCK040	Red medium grit cup composite polisher	30	\$59
FPCCK010	Gray fine grit point composite polisher	15	\$59
FPCCK040	Gray fine grit cup composite polisher	15	\$59



Intra-oral Adjustment/Polishing
for Zirconia and Lithium Disilicate

Adjust and Polish Zirconia and Lithium Disilicate the eZr Way!

eZr Intra-oral Adjustment:



- **No chatter, chipping or microfracture** – Unlike diamonds bonded directly to a hard metal shaft, eZr's rubber substrate cushions the diamonds to prevent damage to the restoration.
- **No heat build-up** – The carefully engineered substrate does not transfer heat during use.
- **FAST!** – eZr's coarse grit blue grinders make your adjustments fast and convenient.

eZr Intra-oral Polishing:



- **Unparalleled smooth polish** – Don't take our word for it! – "The eZr system is the only one I've tried that polishes zirconia to smooth!" Todd Brower, DDS, Holland, Michigan.
- **Increased physical stability** – Highly polished zirconia and lithium disilicate is stable and non-abrasive to opposing dentition.
- **No paste required** – Forget the mess and expense of polishing pastes. eZr™ is all you need.
- **No water required** – No need to use water when polishing. Water is only required for high speed.

eZr™ is a three-step diamond polishing system for the processing of zirconium oxide and lithium disilicate. The grinders and polishers are ideal for adjusting and correcting chairside restorations such as inlays, onlays, veneers and crowns.



*Slow speed right
angle latch shank
Steam autoclaveable*



1

Start by gross adjusting and contouring using the blue coarse grit grinders. Use the flat end taper for occlusal/labial/lingual/buccal surfaces and the disc for interproximal surfaces.



2

Use the green medium grit polishers for smoothing and pre-polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



3

Use the orange fine grit polishers for high gloss polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.

Item	Description	Qty.	Price
FPZK01	Complete system kit Gross removal, smoothing, pre-polishing, and high gloss polishing	2 grinders 6 polishers	\$196
FPZC030	 Blue coarse grit disc Gross removal	5/pk	\$157
FPZC010	 Blue coarse grit flat end taper Gross removal	5/pk	\$118
FPZM010 (not in kit)	 Green medium grit small point polisher Smoothing and pre-polishing	5/pk	\$48
FPZM020	 Green medium grit large point polisher Smoothing and pre-polishing	5/pk	\$74
FPZM030	 Green medium grit disc polisher Smoothing and pre-polishing	5/pk	\$74
FPZM040	 Green medium grit cup polisher Smoothing and pre-polishing	5/pk	\$74
FPZF010 (not in kit)	 Orange fine grit small point polisher High gloss polishing	5/pk	\$48
FPZF020	 Orange fine grit large point polisher High gloss polishing	5/pk	\$74
FPZF030	 Orange fine grit disc polisher High gloss polishing	5/pk	\$74
FPZF040	 Orange fine grit cup polisher High gloss polishing	5/pk	\$74

Cases using other matrix systems



Images courtesy of Dr. Julio Maya



Images courtesy of Dr. Dory Stutman



Images courtesy of Dr. Matthew A. Nejad, DDS, Beverly Hills, CA



Images courtesy of Dr. Jack D. Griffin, Jr., DMD, MAGD, Weldon Spring, MO



Images courtesy of Dr. Jack D. Griffin, Jr., DMD, MAGD, Weldon Spring, MO

Reel Matrix™

TESTIMONIALS

"Using a matrix system like Reel Matrix can enable them (doctors) to quickly and precisely place anatomically correct Class II posterior composite restorations that demonstrate exceptional contacts."

Dr. Dory Stutman

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Dr. Matthew A. Nejad, DDS
Beverly Hills, CA*

"When there is less than 2mm of separation between the adjacent preparations, back-to-back restoring should be considered. In this case, Pinch™ matrices were placed on the first bicuspid and molar and tightened. These feature a blue mylar and an attached tension retainer that is pinched, securing it to the tooth. Slick Band™ sectional matrices are then placed on adjacent restorations and held in place with wedges."


*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*

"To ensure a tight contact, a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument. This ensures that the contact is tight and well-formed while light curing though it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*

PerForm™

Cases using anteriors/cervicals



Images courtesy of Dr. Zunzarren, DDS



Images courtesy of Dr. Abdi Sameni, DDS


VariStrip™

How can you produce a natural curvature for anterior composites and eliminate black triangles? Check out Garrison's curved "Varistrips" to recreate the natural interproximal anatomy for anterior teeth.

- Helps avoid incorrect flat embrasures, black triangles and food traps that can be caused by old-fashioned straight mylars.
- No need for wedges on subgingival fillings, as the contours naturally close off the apical aspects for you.
- Excellent for class IVs, diastema closures, etc.
- Variable height allows you to slide each strip in until it precisely matches the tooth height before cutting off excess.
- The blue tint does not affect cure, but ends the annoyance of losing clear strips on patient bibs or bracket tables.

Blue View VariStrips—a must-have super-simple product for every dental office.

Dr. Michael Curtis, Bridgeport, CT, USA



Clinical Tip

After using your preferred bonding procedure, place the restorative material into the preparation in advance of using the cervical matrix.

Place the cervical edge of the matrix against the tooth (Fig. 1).

While pressing the gingival edge of the matrix firmly against the tooth, roll the matrix onto the tooth (Fig. 2). Excess material will be expressed from the occlusal edge of the matrix rather than subgingivally, greatly simplifying clean up.

Cure directly through the matrix.

ReelMatrix

Real Easy. Real Tight. Real Fast!

For all posterior restorations!

1 Real Easy:
The reels pop into the handle without any fuss! Press the reel into the handle and retract the grippers to lock the reel in place. Plus, if you use the regular matrices (yellow hubs) they can be easily pushed through contacts.

2 Real Tight:
Tight contacts and excellent contour are what Reel Matrix is all about. These thin matrices are pre-contoured in three dimensions and when combined with anatomical wedges and separator rings, the result is amazingly tight contacts at the natural height of contour.

3 Real Fast:
Pre-loaded reels mean you can just lock them in and go. The retainerless design greatly speeds up quadrant work while improving visibility and patient comfort. Add an FX500 separator ring and you virtually eliminate buccal and lingual flash.



Margin Elevation Technique restores deep caries without the trauma.

Margin elevation requires a matrix that will properly adapt and isolate the subgingival area. This is not possible with traditional matrix materials (Figure A). The narrow Margin Elevation Band, with its enhanced conical design, is able to slide down subgingivally and properly seal this area (Figure B).




Figure A




Figure B

For the complete margin elevation technique, go to garrissdental.com



RM500H
Margin elevation



RM100D
5mm deadsoft



RM200D
6.3mm deadsoft



RM100H
5mm regular



RM200H
6.3mm regular



RM100B
5mm blue



RM200B
6.3mm blue

The Reel Matrix™ is the ideal delivery system for the Margin Elevation Band. Its extremely small size and light weight prevents the accidental dislodgement that can occur with long, cumbersome metal retainers.



Margin elevation band is .002" thin.
Super-thin dead-soft bands are .0012" thin.
Regular bands are .0015" thin and can be placed through existing contacts. A Blue View cure-through option at .0029" thin is available if you prefer trans-enamel curing.

ReelMatrix

Real Easy. Real Tight. Real Fast!

Item	Description	Price
RMK05	Kit Contains 80 matrices (20 RM500H, 10 ea. RM100D, RM200D, RM100H, RM200H, RM100B, RM200B) 1 tensioning instrument 1 PerForm™ instruments (set of 2)	\$471

Refills		Qty.	Price
RM500H	Margin elevation band	50	\$91
RM100D	5mm dead-soft matrices	50	\$83
RM200D	6.3mm dead-soft matrices	50	\$83
RM100H	5mm regular matrices	50	\$83
RM200H	6.3mm regular matrices	50	\$83
RM100B	5mm blue matrices	50	\$93
RM200B	6.3mm blue matrices	50	\$93
RMR100	5.9mm empty reel	100	\$53
RMR200	6.7mm empty reel	100	\$53
RMH-100	Tensioning instrument		\$240

IMPORTANT ACCESSORIES

PerForm™ Proximal Contact Instrument
(included in RMK05 kit)

Page 35

Great for core build-ups

SlickBands™

Tofflemire-style

Non-stick technology developed by Garrison®: Modern bonding agents can adhere to uncoated matrix bands—Slick Bands™ slide right out.

Reaches over to the adjacent tooth more easily to help you create broader, deeper contacts.

- Extremely wide preps are simplified.
- Improved visibility for core buildups.



The increased curvature of the all-new **Right-Curve matrix band** produces outstanding gingival seal and reaches over to adjacent teeth easier for broader, deeper contacts. This creates a significantly increased funnel-shaped cone when wrapped around a tooth. It's the right choice for many restorations.

Tighter gingival seal is beneficial for all restorations.

Traditional band with narrow 'funnel'.



Margin Elevation Band


The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.



A margin elevation matrix will give you the best gingival adaptation and contour for the deep margin areas. It is trimmed in height to allow improved adaptation to deep margin areas.

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Matthew A. Nejad, DDS
Beverly Hills, CA*



Slick Bands™ Tofflemire-style Matrices


Clinical Tip: Use dead-soft green bands when contact is COMPLETELY broken with adjacent teeth for maximum burnishability. When contact is NOT broken with adjacent teeth, use regular gray bands.

IMPROVED Dead-soft Bands

with non-stick coating



#1 Universal for general restorations



Right-Curve™ matrices provide a lower contact point and improved isolation. Excellent for all restorations.




#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations

Tofflemire-style matrices are .0014" thick
Margin elevation matrices are .002" thick




Regular Bands


with non-stick coating




#1 Universal for general restorations



Right-Curve™ matrices provide a lower contact point and improved isolation. Excellent for all restorations.



#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations



SMTK06

Slick Bands™

Description
Slick Bands Kit
300 Assorted Slick Bands
100 G-Wedges
PerForm instruments (1 set; shown on page 35)

Item
SMTK06
Kit
\$339

		Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Dead Soft Bands	#1 Universal Tofflemire	50	\$41	100	\$65
	Right-Curve	50	\$52	100	\$80
	#13 Pedodontic	50	\$41	100	\$65
	#2 Sub-gingival	50	\$41	100	\$65

		Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Regular Bands	#1 Universal Tofflemire	50	\$37	100	\$58
	Right-Curve	50	\$49	100	\$77
	#13 Pedodontic	50	\$37	100	\$58
	#2 Sub-gingival	50	\$37	100	\$58

		Standard (-M)		Bulk	
		Qty.	Price	Qty.	Price
Regular Bands	Margin elevation band	50	\$57	100	\$85

To order standard packs, add a -M after the item #.

Pinch[™] Matrices

Your fastest and easiest choice when you need a Tofflemire-style matrix.

Not every Class II can be restored with Composi-Tight®, so when you need a full matrix, Garrison® has you covered. Pinch[™] Matrices are retainerless Tofflemire-style matrices, in cure through or metal, with a built-in tensioning ring. There's no more fiddling with retainers or fancy tightening devices, simply slip around the tooth to be restored and Pinch!

Retainerless:



The easiest-to-place circumferential matrix eliminates bulky retainers.

Pre-contoured:

Anatomy-hugging shape helps produce properly contoured restorations.

Blue Tint:


Provides contrast between matrix and tooth structure without compromising composite polymerization.



Pinch Matrices

Metal Kit	Description	Qty.	Price
KMMK1	Metal Assortment Kit 50/ea. size	150	\$233
Metal Refills			
KMM100	Metal Premolar	50	\$79
KMM175	Metal Molar	50	\$79
KMM200	Metal Tall Molar	50	\$79
Blue View Kit			
KPMK1	Cure Through Assortment Kit 50/ea. size	100	\$168
Blue View Refills			
KPM100	Cure Through Premolar	50	\$85
KPM200	Cure Through Molar	50	\$85


How PerForm™ works:



1A



1B



2

STEP 1:

Place the first composite increment. Press PerForm™ into the uncured composite increment and either push (distal; fig. 1A) or pull (mesial; fig. 1B).

STEP 2: Place the curing light directly on the Fresnel lens and cure.

PerForm

Proximal Contact Instruments

PerForm puts you in control.

- Improves contacts with any matrix system
- Control the height and breadth of contacts
- Increase the tightness of contacts
- Unique lens controls curing light diffusion
- Steam autoclave only

PerForm™

Item	Kit	Qty.	Price
CFK01	PerForm (1/ea. size)	2/pk	\$122
CF100	1.1 mm Small Instruments	2/pk	\$122
CF200	2 mm Large Instruments	2/pk	\$122



"To ensure a tight contact a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument. This ensures that the contact is tight and well-formed while light curing through it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*



Photo
courtesy of
David S. Hornbrook, DDS, FAACD
International Lecturer and
Private Practice, San Diego, CA


Curved Form/tapered: Provides the optimum curvature and band height for virtually every anterior restoration.

Blue Tint: Provides contrast between matrix and tooth structure without compromising composite polymerization.

Our Most Popular Anterior Matrix

BlueView™ VariStrip™

Contoured Anterior Matrices



The Blue View™ VariStrip™ is a revolutionary contoured anterior matrix providing the optimum curvature and band height for virtually every anterior restoration. The 0.0020" thin plastic anatomical strip is tapered from one end to the other. You can place the strip interproximally and then slide it until the tooth height is perfectly matched. The pre-contouring easily recreates the occluso-gingival anatomy, avoiding flat embrasures.

Blue View™ VariStrip™

Item	Kit	Qty.	Price
VS01	VariStrip	100	\$115

BlueView™

Cervical Matrices

"I love the Blue View Cervical Matrix! It works so well. My assistants love watching me use this system. The patients love the excitement in our faces. It makes the process easy and fast. Gotta love that! We are also hooked on using the ReelMatrix. Thanks so much for making great products."

Cynthia Brattesani, DDS,
San Francisco, CA, USA



Flexible

- adapts to individual anatomy for improved contours
- use it to directly compact the composite into the restoration
- thin edges allow excess composite to flow out for improved anatomy

Blue Tint


- provides contrast between matrix and tooth structure
- does not compromise composite polymerization

Convenient!

- protects the restoration from contamination
- eliminates time spent hand sculpting
- eliminates air-inhibited layer
- reduces finishing time
- 7 sizes for a perfect fit
- Gingival Retractor included!

KCMK1 INCLUDES








- 275 Assorted Matrices
- 1 Matrix placement Instrument
- 1 Gingival Retractor



MATRICES INCLUDED:

- Anterior/Premolar
KCMA60 - 50 pcs
- Molar
KCMA70 - 50 pcs
- KCMA10 - 35 pcs
- KCMA20 - 35 pcs
- KCMA30 - 35 pcs
- KCMA40 - 35 pcs
- KCMA50 - 35 pcs

Blue View™ Cervical Matrices

Item	Description	Price
KCMK1	Assortment Kit w/Instruments	\$184
Refills		
KCMA10	 — Refill 100 pack	\$102
KCMA20	Refill 100 pack — 	\$102
KCMA30	 — Refill 100 pack	\$102
KCMA40	Refill 100 pack — 	\$102
KCMA50	 — Refill 100 pack	\$102
KCMA60	Refill 100 pack — 	\$102
KCMA70	 — Refill 100 pack	\$102
KCMINST1	Placement Instrument	\$21

Composi-Tight® Matrix Forceps


Precise Placement, Unbeatable Grip.



Check out the GRIP! Effortlessly remove the matrix bands that otherwise would be a struggle.



Precise Placement: The tips have a 75-degree angle with an occlusal-gingival orientation, providing a more appropriate placement angle than cotton forceps. The petite, carbide tips allow the matrix to be grasped at its upper edge to prevent distortion of the interproximal area of the band.



Unbeatable Grip: The tungsten carbide tips provide exceptional grip and tactile control for tough access cases.

Matrix Forceps

BFA-1 Angled tip for occlusal placement \$155



Precision Instrumentation completes the NiTin™ system!

Low profile tips simplify stacking on MODs with NiTin's smaller rings. Recommended for NiTin™ rings only.



This surgical tissue forceps with its deep interdigitating teeth provide ample grip and tactile control for both band and wedge placement.

NiTin™ Precision Instruments

NTG-MRDF-100	Matrix Forceps	\$77
NTG-DF6-8	Surgical Tissue Forceps	\$26

toll free 888.437.0032 | www.garrisondental.com

Composi-Tight® 3D Fusion™

Composi-Tight® 3D Fusion™ includes dentistry's premier heavy-duty stainless steel ring placement forceps. Regardless of your hand size or strength, Composi-Tight® 3D Fusion™ forceps will give you incredible control for the precise placement of super-strong separator rings over even the widest molars.


**Composi-Tight® 3D Fusion™
Ring Placement Forceps**
FXP01 \$99

Composi-Tight® 3D Fusion™ forceps will work with virtually any ring from any manufacturer.



Ring Placement Forceps

Holds all your sectional matrix system rings.



Ring Placement Forceps

AUMRDF-100 \$83

Universal Forceps

for placing rubber dam
clamps and rings

UNVF-100 \$83



Composite Instruments Nonstick titanium-nitride coated

Our multi-function instruments have double condenser configurations and carving blades on opposite ends. Now you can place, pack, shape and carve your way to perfectly contoured composite restorations without ever having to switch instruments!

- Large-diameter, lightweight handles for improved ergonomics.
- Provides non-stick surface for placement and shaping of composites.
- 40% smoother than stainless steel for a non-porous, non-stick surface.
- Withstands diverse cleaning/sterilization methods, including ultrasonic cleaning.



“Skippy Special”
TN008

Item	Description	Price
TN008	“Skippy Special” 2 condensers, ball burnisher, blade, acorn	\$105



5-in-1
TN009

Item	Description	Price
TN009	5-in-1 Multi-Function Instrument 2 blades, 2 condensers, acorn	\$105

Item	Description	Price
TN001	Medium Blades Placement of composite material Gingival retraction to sculpt under free margins Burnishing matrix bands	\$45



TN002	Medium Blade with Medium Plugger Contouring small Class I, II, V Packing composite material into conservative preparations Burnishing matrix bands	\$45
-------	---	------



TN003	Medium Blade with Large Plugger Contouring large Class I, II, V Packing composite material into large preparations Convex plugger tip minimizes pullback of composite material	\$45
-------	---	------



Item	Description	Price
TN004	Extra Thin Blades Reach under margins to hone and sculpt Sculpting interproximally Packing cord in gingival crevice	\$48



TN005	Long Extra Thin Blades Allows additional access for hard-to-reach posterior contouring and shaping Sculpting interproximally Packing cord in gingival crevice	\$48
-------	--	------








TN006	Small/Medium Multi-Function Plugger Placement of composite material in smaller Class II and III Packing material toward gingival seats Carving occlusal anatomy with acorn carver	\$73
-------	--	------



Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials	Page	Special Notes
Distal of canine, short or malpositioned 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! All new separator rings designed to hold firmly on tough cases.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Conservative Class II 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Slightly wide and/or deep 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.

Procedure	Recommended Materials	Page	Special Notes
Large Class II with missing cusps 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit—Green Wide Prep Ring	4–5	EXCLUSIVE! Dentistry's first Wide Prep ring greatly simplifies these more challenging restorations.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
Small to medium size Class III 	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Large Class III or Class IV 	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Class V 	Blue View™ Cervical Matrix	37	Contours and directly compresses resin while eliminating air-inhibited layer.
Crown build-ups, post and core 	Reel Matrix™	30–31	Best choice. Margin Elevation band available for very deep buildups. Retainerless convenience.
	Pinch™ Matrices—Metal or Blue View	34	Very fast; may not provide as tight a cervical seal as other matrices.
	Slick Bands™ Tofflemire-Style Non-Stick Matrices	32–33	Margin Elevation band available for very deep buildups if you prefer a traditional matrix retainer.

Order Form

Fax: 616.842.2430

Item Number	Description	Quantity	Price/Unit	Total Price
		Shipping		
		Sales Tax		
		Total Price		

Call for additional express shipping rates and delivery. All prices subject to change without notice. Appropriate shipping charges and sales tax will be added. Special offers not valid with any other advertised specials.

Doctor's Name

Contact Name

Address

City State Zip

Phone _____ Fax _____

Email

Payment Method

Credit Card: ☐ Visa ☐ MasterCard ☐ Am. Express ☐ Discover ☐ C.O.D. (add \$13.00)

Credit Card #	CC Ref #
---------------	----------

Exp. Date	Signature
-----------	-----------

Credit Card Billing Information

Address

City State Zip

toll free 888.437.0032 | garrisondental.com

Customer Service Information

Call **888.437.0032 Toll Free**
(USA & Canada)

616.842.2244

Fax 616.842.2430

Email gds@garrisondental.com

Online **garrisondental.com**
Look for our Learning Center
with clinical tips, case studies
and technique guides.

Mail 150 DeWitt Lane
Spring Lake, MI 49456
USA

Shipping Rates

Regular UPS	2nd Day Air	Next Day Air	C.O.D.
\$9.75	\$19.00	\$34.00	+\$13.00

Express orders must be placed by 4pm EST to ship same day.

Shipping charges subject to change without notice.

Payment Methods

VISA, MasterCard, American Express,
Discover, Check Prepay or C.O.D.

Prices are subject to change without notice.

Distribution

Garrison is a global leader for practical and innovative dental products and accessories, with distributors worldwide.

We are pleased to announce that your favorite Garrison products are available through several full-service dealers in the USA.

To find a dealer, contact our headquarters in Spring Lake, Michigan. Phone: 616.842.2244 or Email: gds@garrisondental.com



Order Toll-Free • Order by Fax • Order Online • Order by Mail



Our Guarantee – It's Simple!

Satisfaction = 100%

Everything you order from Garrison® Dental Solutions is backed by our risk-free six-month money-back **Satisfaction Guarantee**. If you are not satisfied with your product, return the unused portion and we will issue a refund.

Quality = 100%

All Garrison products are covered by our **Quality Guarantee**. We will repair or replace any Garrison® product that is shown to be defective in materials or workmanship based on the individual product's anticipated lifetime. Please visit our website for details on our various products, garrisondental.com.

toll free 888.437.0032 | www.garrisondental.com

INCAT121