

2021 PRODUCT CATALOG



# table of contents

## **MATRIX SYSTEMS**

sectional matrix systems	Composi-Tight® 3D Fusion™ . . . . .	4-5
	Other sectional systems and refills . . . . .	8-9
retainerless Tofflemire-style	Reel Matrix™ . . . . .	30-31
	Pinch™ Matrix . . . . .	34
traditional Tofflemire-style	Slick Bands™ . . . . .	32-33
anterior matrix systems	VariStrip™ . . . . .	36
cervical matrix systems	Blue View™ Cervical Matrices. . . . .	37

## **WEDGES** plastic

Composi-Tight® 3D Fusion™ Wedge . . . . .	14
A+Wedge® . . . . .	15
WedgeWands® . . . . .	15
G-Wedge™ . . . . .	16

## wood

SoftWedge™ . . . . .	16
----------------------	----

## pre-wedge & shield

FenderWedge® . . . . .	17
------------------------	----

## **IMPORTANT ACCESSORIES**

polishing systems	Rally™ . . . . .	24-25
	eZr™ . . . . .	26-27
interproximal finishers	FitStrip™ . . . . .	20-23
forceps	Composi-Tight® 3D Fusion™ Forceps . . . . .	39
	Composi-Tight® Matrix Forceps . . . . .	38
	NiTin™ Forceps . . . . .	38
	Ring Placement Forceps . . . . .	39
	Universal Placement Forceps . . . . .	39
specialty	PerForm™ Contact Former . . . . .	35
composite instruments	Multi-Function Non-Stick . . . . .	40
	Standard Non-Stick . . . . .	41


## **INFORMATION PAGES**

Composi-Tight® 3D Fusion™ . . . . .	2-3, 6-7
Band Comparison Chart . . . . .	10-12
Matrix Selection Guide . . . . .	42-43
Order Form . . . . .	44
Ordering Information . . . . .	45
Guarantee . . . . .	back cover

**Only a sectional matrix can deliver results like this!**

### Tofflemire-style System

- ❌ Fails to restore proximal anatomy
- ❌ Thin contact at the marginal ridge
- ❌ Large food trap below
- ❌ Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



### Composit-Tight® 3D Fusion™ Sectional Matrix System

- ✓ Operator-friendly retaining system
- ✓ Naturally contoured bands
- ✓ Anatomically correct contacts
- ✓ Contacts at the height of contour
- ✓ Rings produce optimal tooth separation for excellent, tight contacts

**Only Composit-Tight® 3D Fusion™ can make the technique this intuitive and frustration free!**



#### 1 Place the matrix band:

Dead-soft three dimensionally contoured matrices have placement tabs for easier occlusal placement and Garrison's revolutionary non-stick coating for hassle-free removal.



#### 2 Wedge firmly:

Composit-Tight® 3D Fusion™ Ultra-Adaptive wedges glide in but won't back out. Their soft, flexible exterior readily adapts to interproximal irregularities for a superior gingival seal.



#### 3 Place the ring:

All three Composit-Tight® 3D Fusion™ separator rings will fit right over the wedge for easier placement. With Ultra-Grip™ retention extensions, they will stay firmly in place from the distal of the canine to the widest of molar preparations.

# Composi-Tight® 3D Fusion™ Cases

Great adaptation to the teeth eliminates the need for adjustment with polishers in the proximal areas.



Images courtesy of Clarence Tam, Auckland, NZ

Wide preps can be handled in one simple step when you have the right (FX600 green) ring.



Images courtesy of Dr. Nikhil Sethi, BDS (Lon) Msc (Lon), Birmingham, UK

Back-to-back requires a shift in technique—we can show you how.\*



Images courtesy of Dr. Jeffrey Kelley, DDS, Zeeland, MI

Quadrants can be a challenge. Ease your concerns with tips and tricks from the pros.



Images courtesy of Dr. Mike Malone, DDS, Lafayette, LA

\*To see full writeups, cases and reviews go to [www.garrisondental.com/learning-center](http://www.garrisondental.com/learning-center)  
Do you have a case you want to share? Email it to us at [contact@garrisondental.com](mailto:contact@garrisondental.com).

## The Award Winning Composi-Tight® 3D Fusion™ Sectional Matrix System



### TESTIMONIAL

"3D rings are the real topic of Garrison's systems. The "v" shape of the ring that fits in the interproximal area allows a good fitting between the cavity margins and the matrix in the buccal and palatal walls. This results in an easier positioning of the composite masses close to the cavity margins and also the final remodeling usually necessary at the time of removal of the matrix will be very minimal. The rings also permit a divergence of the interproximal dental elements, and this causes a great point of contact. Garrison systems make Class II restorations simpler and more predictable and also reduce the operating time of the finishes when the matrix is taken off."

Enrico Cogo, DDS  
Legnago, Ferrara, and Goito, Italy

# Composi-Tight® 3D Fusion™



**A system designed  
to handle even the  
toughest cases**

More than 20 years of sectional matrix leadership and an incredible amount of feedback from users all over the globe were poured into the design. All this hard work has resulted in separator rings with significantly improved performance and ease of use.



**Hugs the curves** – Soft-Face™ silicone tips with enhanced marginal ridge anatomy hug the matrix band to the tooth for flash free restorations

**It just keeps going...** Enhanced durability for extended ring life and greater tooth separation

## Wider indications, wider smiles!



Dentistry's first Wide Preparation separator ring (green) makes Composi-Tight® 3D Fusion™ the most versatile and user-friendly sectional matrix system available. The Wide Prep ring greatly simplifies what until now has been one of the most challenging posterior composite restorations. Now you can restore these big preps with confidence knowing that you will achieve excellent contact and contour without having to take extra steps.

## Matrix band perfection

Composi-Tight® 3D Fusion™ Full Curve Non-Stick Matrices are three dimensionally contoured, ultra-thin and dead-soft to help you recreate ideal interproximal anatomy.



**Easy in – easy out!** A Grab-Tab™ centered on the occlusal edge of the band simplifies placement while Slick Bands™ non-stick coating reduces composite adhesion by an impressive 92%!

**Go deep or wide** – Improved subgingival extensions and added wrap-around length make for the perfect matrix band to use with the green Wide Preparation ring. Now you're ready to tackle even the widest preps!



FXR-KFF-11

See a product demo online!



# Composi-Tight® 3D Fusion™

Sectional Matrix System



Composi-Tight® 3D Fusion™ Bands # .0016" thick	FX100	FX150	FX175	FX200	FX300
Non-stick dead soft stainless steel	bicuspid 4.0 mm w/ext 4.4 mm	bicuspid with deep cervical preparations 4.1 mm w/ext 6.0 mm	bicuspid and small molars 5.0 mm w/ext 5.6 mm	molars 5.4 mm w/ext 6.6 mm	molars with deep cervical preparations 6.0 mm w/ext 8.7 mm

### Kit

- 1 Short Ring (blue)
- 1 Tall Ring (orange)
- 1 Wide Prep Ring (green)
- 70 Assorted Full Curve Non-Stick Matrix Bands
- 80 Assorted Ultra-Adaptive Wedges
- 4 Assorted FenderWedges® (pg 17)
- 12 Assorted Rally™ Polishers (pg 24–25)
- 1 Ring Placement Forceps (pg 39)
- 1 BFA-1 Band Placement Forceps (pg 38)
- 1 TN009 Multi-Function Instrument (pg 40)

### Item

- FXR-KFF-11
- FXR-KFF-10 (BFA-1 and TN009 not in kit)
- FXR-KFF-14 (TN009 not in kit)

### Composi-Tight® 3D Fusion™ Rings

- Short Ring (blue)
- Tall Ring (orange)
- Wide Prep Ring (green)
- 3-Ring Kit (1 of each ring)

Item	Qty.
FX400	2
FX500	2
FX600	2
FXR01	3

### Composi-Tight® 3D Fusion™ Bands

- Bicuspid
- Bicuspid w/ deep cervical preps
- Bicuspid and small molars
- Molars
- Molars w/ deep cervical preps

Item	Standard (-M)	Bulk
Qty.	Qty.	Qty.
FX100	50	100
FX150	30	60
FX175	50	100
FX200	50	100
FX300	30	60



## Case Study

*"The Composi-Tight® 3D Fusion™ system allows you to achieve excellent results on any clinical condition. The Composi-Tight® 3D Fusion™ Firm bands, in particular, allow you to gently move the rubber dam and tissues cervically, resisting deformation during positioning. The combined use with Composi-Tight® 3D Fusion™ rings allows you to stabilize and adapt the matrix to the anatomy of the teeth and to recreate a tight contact point."*

Dr. Elio Boschetti, Arzignano (VI)  
ESCD Certified Member



**Fig. 1** Preoperative view with subtle dark caries visible through the translucent enamel and cavities prepared and finished.



**Fig. 2** Sectional matrix system "ring-band-wedge" in order to achieve a natural anatomy of inter-proximal walls and correct contact areas.



**Fig. 3** Final posterior composite resin restorations after finishing and after high gloss polishing.



## Case Study

*"This is an ideal case for Garrison® Composi-Tight® 3D Fusion™ Green FX600 separator rings: Distal Class II on tooth 15. The different profile and special geometry of the tips allows an ideal separation and a correct profile at interproximal level between teeth with different morphologies."*

Dott. Andrea Fabianelli  
EPA European Prosthodontist Specialist  
Studio dentistico associato Fabianelli Plahuta



**Fig. 1** Class II lesion , tooth 15, detectable with transillumination or x-ray.



**Fig. 2** Opening of cavity, which appears quite extensive.



**Fig. 3** Cavity as it appears at the end of the removal of infected tissues and with well defined margins.



**Fig. 4** Positioning of sectional matrix and separator ring: the different morphology of its tips allows a correct modeling of the interproximal wall.



**Fig. 5** Final view of restoration.



**Fig. 6** Restoration after removal of the rubber dam and occlusal check.

# Composi-Tight® 3D XR

Sectional Matrix System



Slick Bands™ XR .0016" thick	SXR100	SXR150	SXR175	SXR200	SXR300
Non-stick dead soft stainless steel	bicuspids 4.6 mm	bicuspids with deep cervical preparations 3.8 mm w/ext 5.4 mm	bicuspids and small molars 5.4 mm	molars 6.4 mm	molars with deep cervical preparations 6.4 mm w/ext 8.6 mm

## Kit

### 3D XR System Kit

- 1 Composi-Tight® 3D XR Ring
- 2 Original 3D500 Soft-Face™ Rings
- 100 Slick Bands™ XR Non-stick Matrices  
(30 ea. SXR100, SXR200, SXR175; 5 ea. SXR150, SXR300)
- 100 WedgeWands® (25 ea. size)
- 20 FenderWedges® (5 ea. size)
- 1 Ring Placement Forceps (AUMRDF-100))

## Item

3D-KSW-10

## 3D Rings

Soft-Face™ Extra Retention Rings  
Soft-Face™ Rings

Item	Qty.
3DXR	2
3D500	2

Composi-Tight® SXR Bands	Item	Standard (-M)	Bulk
		Qty.	Qty.
Bicuspids	SXR100	50	100
Bicuspids w/ deep cervical preps	SXR150	30	60
Bicuspids and small molars	SXR175	50	100
Molars	SXR200	50	100
Molars w/ deep cervical preps	SXR300	30	60

# Composi-Tight® 3D Clear

Cure-through sectional matrix system



Composi-Tight® 3D Clear .0020" thick	MC100	MC175	MC200
Blue tinted polyester	premolar 5.0 mm	bicuspids and small molars 5.9 mm	tall molar 6.8 mm

## Kit

### 3D CLEAR System Kit

- 1 Composi-Tight® 3D XR Ring
- 2 3D600 Clear Soft-Face™ Rings
- 100 Cure-Through matrices  
(30 ea. MC100, MC200; 40 ea. MC175)
- 100 Cure-Through WedgeWands® (25 ea. size)
- 20 FenderWedges® (5 ea. size)
- 1 Ring Placement Forceps (AUMRDF-100))

## Item

3D-KCC-10

## 3D Clear Rings

Soft-Face™ Extra Retention Rings  
Clear Soft-Face™ Rings

Item	Qty.
3DXR	2
3D600	2

3D Clear Bands	Item	Standard (-M)	Bulk
		Qty.	Qty.
Premolar matrix bands	MC100	50	100
Bicuspid and small molar matrix bands	MC175	50	100
Tall molar matrix bands	MC200	50	100

To order standard packs, add a -M after the item #.  
Complete band selection guide pages 10, 11 and 12.

# Composi-TightGold®

## Sectional Matrix System



Composi-Tight Gold # .0015" thick	AU050	AU100	AU150	AU200	AU300
Gold colored dead soft stainless steel	pediatric 3.2 mm w/ext 4.2 mm	bicuspid 4.6 mm	bicuspid with deep cervical preps 3.8 mm w/ext 5.4 mm	molars 6.4 mm	molars with deep cervical preps 6.4 mm w/ext. 8.6 mm

### Kit

#### Gold System w/Forceps

- 250 Matrix Bands  
(100 ea. AU100, AU200; 25 ea. AU150, AU300)
- 6 G-Rings (3 standard, 3 long)
- 1 Ring Placement Forceps

### Item

AUK2

### Gold G-Rings™

G-Rings™ w/Standard Length Tines

G-Rings™ w/Long Length Tines

### Item

AU400 3  
AU500 3

Gold Bands	Item	Standard (-M)	Bulk
		Qty.	Qty.
Pediatric	AU050	25	50
Bicuspid	AU100	50	100
Bicuspid w/ deep cervical preps	AU150	25	50
Molars	AU200	50	100
Molars w/ deep cervical preps	AU300	25	50

### Forceps

Ring Placement Forceps



### Item

AUMRDF-100

# Composi-Tight™

## Sectional Matrix System



Composi-Tight # .0015" thick	B050	B100	B150	B200	B300
Polished dead soft stainless steel	pediatric 3.2 mm w/ext 4.2 mm	bicuspid 4.6 mm	bicuspid with deep cervical preps 3.8 mm w/ext 5.4 mm	molars 6.4 mm	molars with deep cervical preps 6.4 mm w/ext. 8.6 mm

### Kit

#### Original System w/Forceps

- 225 Matrix Bands (100 ea. B100, B200; 25 B300)
- 7 G-Rings (4 standard, 3 long)
- 1 Ring Placement Forceps

### Item

K102

### Original G-Rings™

G-Rings™ w/Standard Length Tines

G-Rings™ w/Long Length Tines

### Item

G100 3  
G200 3

Original Bands	Item	Standard (-M)	Bulk
		Qty.	Qty.
Pediatric	B050	25	50
Bicuspid	B100	50	100
Bicuspid w/ deep cervical preps	B150	25	50
Molars	B200	50	100
Molars w/ deep cervical preps	B300	25	50

To order standard packs, add a -M after the item #.  
Complete band selection guide pages 10, 11 and 12.

# Composi-Tight® 3D Fusion™

**NEW!**

## Firm Matrix Bands

New **Composi-Tight® 3D Fusion™ Firm** matrix bands with **eZ-Place™** technology maintain their shape during placement making them *significantly easier* to insert!

eZ-Place™ technology combines specially processed stainless steel with a proprietary matrix band forming process to yield Fusion Firm's remarkable handling characteristics. They resist deformation nearly 70% better than typical dead-soft bands yet are pre-contoured and super thin at 0.0015". You're ensured the tight contacts and ideal anatomy you've come to expect from the Composi-Tight® 3D Fusion™ system.

### Fusion Firm:

- Nearly 70% better at resisting deformation!
- Super-thin at only 0.0015"
- Five ideal sizes!

### Perfect for:

- Tight interproximal spaces!
- More conservative preps!
- First-time users and veterans alike!



FXH Series

Composi-Tight®  
3D Fusion™ Firm

thickness  
0.0015"

height

Pediatric

Bicuspid

Bicuspid  
with deep  
cervical preps

Bicuspid  
and small molars

Molars

Molars  
with deep  
cervical preps

NA

FXH100  
4.3 mm

FXH150  
4.3 mm  
w/ext 5.9 mm

FXH175  
5.3 mm


FXH200  
6.1 mm

FXH300  
6.1 mm  
w/ext 7.7 mm

Description	Item	Qty.
Assorted Band Kit: 5 sizes	FXHB04	300
Assorted Band Mini Kit: 5 sizes	FXHB05	150

Item	Standard (-M)	Bulk
	Qty.	Qty.
FXH100	50	100
FXH150	50	100
FXH175	50	100
FXH200	50	100
FXH300	50	100

FX Series  
SXR Series  
B Series  
M – SM – AU Series



## Band Comparison Chart

FX Series

Composi-Tight®  
3D Fusion™

Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.

thickness  
0.0016"

NA



FX100  
4.0 mm  
w/ext 4.4 mm



FX150  
4.1 mm  
w/ext 6.0 mm



FX175  
5.0 mm  
w/ext 5.6 mm



FX200  
5.4 mm  
w/ext 6.6 mm



FX300  
6.0 mm  
w/ext 8.7 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
FX100	50	100
FX150	30	60
FX175	50	100
FX200	50	100
FX300	30	60

SXR Series

Slick Bands™ XR

Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.

thickness  
0.0016"

height

NA



SXR100  
4.6 mm



SXR150  
3.8 mm  
w/ext 5.4 mm



SXR175  
5.4 mm



SXR200  
6.4 mm



SXR300  
6.4 mm  
w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
SXR100	50	100
SXR150	30	60
SXR175	50	100
SXR200	50	100
SXR300	30	60

SM Series

Slick Bands™

Slick Bands™ Matrices have non-stick coating.

thickness  
0.0016"

height



SM050  
3.2 mm  
w/ext 4.2 mm



SM100  
4.6 mm



SM150  
3.8 mm  
w/ext 5.4 mm



SM175  
5.5 mm



SM200  
6.4 mm



SM300  
6.4 mm  
w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
SM050		60
SM100	50	100
SM150	30	60
SM175	50	100
SM200	50	100
SM300	30	60







To order standard packs, add a -M after the item #.

## AU Series

	Pediatric	Bicuspsids	Bicuspsids with deep cervical preps	Bicuspsids and small molars	Molars	Molars with deep cervical preps
<b>Composi-Tight Gold®</b>	The Gold bands are the same as the Composi-Tight® bands.					
thickness 0.0015"				NA		
height	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm		AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
AU050	25	50
AU100	50	100
AU150	25	50
AU200	50	100
AU300	25	50

## M Series

<b>Composi-Tight®</b>	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.					
thickness 0.0015"						
height	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm




Item	Standard (-M)	Bulk
	Qty.	Qty.
M050	25	50
M100	50	100
M150	25	50
M175	50	100
M200	50	100
M300	25	50

## B Series

<b>Composi-Tight® Original</b>	The original bands have a slightly shorter length.					
thickness 0.0015"				NA		
height	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm		B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
B050	25	50
B100	50	100
B150	25	50
B200	50	100
B300	25	50

## MC Series

<b>Composi-Tight® 3D Clear</b>	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.					
thickness 0.0020"		NA			NA	
height	MC100 5.0 mm		MC175 5.9 mm	MC200 6.8 mm		

Item	Standard (-M)	Bulk
	Qty.	Qty.
MC100	50	100
MC175	50	100
MC200	50	100

To order standard packs,  
add a -M after the item #.

## Cases using Garrison® Wedges



Images courtesy of Dr. Mike Malone, DDS, Lafayette, LA



Images courtesy of Dr. Robert Rosenfeld, DDS, Westwood, NJ

### Composit-Tight® 3D Fusion™ Wedge

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.



Images courtesy of Dr. Robert Rosenfeld, DDS, Westwood, NJ

### Wedge Wands®

You can twist the handle off or leave it on. If you choose to leave the handle on, it can be used to snug the wedge in tighter once the ring separates the teeth. Then remove it when finished without ever having to pass an instrument back and forth with your assistant.



Images courtesy of Dr. Schuyler Van Gorden, DDS, Eau Clair, WI



Images courtesy of Dr. Strassler, DDS

### FenderWedge®

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations. FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and pre-separates the teeth, thus simplifying the application of a matrix.



Images courtesy of Dr. David S. Hornbrook, DDS, FAACD, San Diego, CA

## TESTIMONIALS

"Garrison's redesigned **3D Fusion™ Wedges** feature soft "fins" that smoothly adapt matrix bands to tooth irregularities and prevent the wedge from backing out. These are helpful in a wide array of situations, to eliminate the likelihood of flash at the gingival margin."

*Dr. Robert Rosenfeld, DDS  
Westwood, NJ*

"The Garrison **A+Wedge™** minimized or eliminated bleeding. It's particularly beneficial on patients with slight to mild gingivitis."

*Joe Arbutante, DDS  
Garrison Advisory Board Member*


"Even those with steady hands and stellar technique occasionally and inadvertently nick the proximal tooth or restoration as they are preparing the adjacent tooth. I know I have! The **FenderWedge®** is a simple device that I wish I would have had the opportunity to use since the day I started practicing dentistry. It would have saved me much aggravation and time. It is now a "must-have" item during all my restorative procedures."

*David S. Hornbrook, DDS, FAACD  
International Lecturer and  
Private Practice, San Diego, CA*

# Composi-Tight® 3D Fusion™

## Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



**A real softie** Composi-Tight® 3D Fusion's Soft-Face™ over-mold allows the wedge to do what no other wedge can truly do – actually adapt to interproximal irregularities.

### Tough on the inside

The firm inner skeleton allows for easy insertion and tooth separation just like a traditional wedge.



### Fins with serious function

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.



## Composi-Tight® 3D Fusion™ Wedge

### Description

Composi-Tight® 3D Fusion™ Wedge Kit  
Composi-Tight® 3D Fusion™ Wedge Kit

### Item

FXK4  
FXK4-M

### Qty.

400  
200

### Refills

X-Small Yellow  
Small Blue  
Medium Orange  
Large Green

### Item

FXYL-M  
FXBL-M  
FXOR-M  
FXGR-M

### Standard Qty.

50  
50  
50  
50

### Item

FXYL  
FXBL  
FXOR  
FXGR

### Bulk Qty.

100  
100  
100  
100

# A+Wedge®

Controls Gingival Bleeding



- Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva
- Improves clarity of restoration area
- Eliminates need to apply an astringent separately

## A+Wedge®

Description	Item	Qty.
A+Wedge® Kit	GWAK4	200

Refills	Item	Qty.	Bulk
X-Small Yellow	GWAYL	100	
Small Blue	GWABL	100	
Medium Orange	GWAOR	100	
Large Green	GWAGR	100	

### Insert



### Wait & Restore



### Remove & Rinse



# WedgeWands®



WedgeWands  
also available in  
Cure-Through

## WedgeWands®

Description	Opaque / Clear	Qty.
WedgeWands® Kit	WK4 / WCK4	400

Refills	Opaque / Clear	Qty.	Standard	Opaque / Clear	Qty.	Bulk
X-Small Yellow	SWYL / SWCYL	100		WYL / WCYL	300	
Small Blue	SWBL / SWCBL	100		WBL / WCBL	300	
Medium Orange	SWOR / SWCOR	100		WOR / WCOR	300	
Large Green	SWGR / SWCGR	100		WGR / WCGR	300	




Perfect wedge. Perfect contact.

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.

## G-Wedge™

Description	Item	Qty.		Item	Qty.
G-Wedge™ Kit	GWK4	400			
Refills	Item	Standard Qty.		Bulk Item	Bulk Qty.
X-Small Yellow	GSWYL	100		GWYL	300
Small Blue	GSWBL	100		GWBL	300
Medium Orange	GSWOR	100		GWOR	300
Large Green	GSWGR	100		GWGR	300



## SoftWedge™

Superior adaptation and seal



**The first** wooden wedge designed for use with a sectional matrix system.

Soft Wedges are wooden interproximal wedges that have been optimized for use with modern restorative materials. The ability to compress during insertion, then swell to help fill voids ensures that your preparations stay clean and dry—an absolute must for today's bonded restorations!

## Soft Wedge™

Description	Item	Qty.		Item	Qty.
Soft Wedge™ Kit	WDK4	400			
Refills	Item	Standard Qty.		Bulk Item	Bulk Qty.
X-Small Yellow	WSDYL	100		WDYL	300
Small Blue	WSDBL	100		WDBL	300
Medium Orange	WSDOR	100		WDOR	300
Large Green	WSDGR	100		WDGR	300

# FENDERWEDGE®

## Interproximal Tooth Shields

### Pre-separate and Protect

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations.

FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and separates the teeth, thus simplifying the application of a matrix.



"Even those with steady hands and stellar technique occasionally and inadvertently nick the proximal tooth or restoration as they are preparing the adjacent tooth. I know I have! The Fender Wedge is a simple device that I wish I would have had the opportunity to use since the day I started practicing dentistry. It would have saved me much aggravation and time. It is now a "must-have" item during all my restorative procedures."

David S. Hornbrook, DDS, FAACD  
International Lecturer and  
Private Practice, San Diego, CA

### FenderWedge®

**Description**  
FenderWedge® Kit

**Item**    **Qty.**  
TSK01    80

Refills	Item	Standard	Item	Bulk
		Qty.		Qty.
X-Small Purple	TSPR-M	30	TSPR	60
Small Orange	TSOR-M	30	TSOR	60
Medium Green	TSGR-M	30	TSGR	60
Large Yellow	TSYL-M	30	TSYL	60



FenderWedge is NOT a sectional matrix system nor does it replace the use of a matrix band.

# Great technique tips and product help are only a click away!

*Join thousands of subscribers on our YouTube channel and get notified instantly when new content is posted.*

The Garrison® YouTube channel is NOT a collection of infomercials. It's all about helpful information on a wide variety of topics. From short two-minute tips on how to reduce flash on Class IIs, to more in depth live patient case presentations, the Garrison® channel delivers great content.

Go to [youtube.com/garrisondental](https://youtube.com/garrisondental) or search Garrison Dental on youtube and subscribe today!



*Speaking of great content...*

## Get CE credits direct from Garrison®!

Go to the Learning Center at [garrisondental.com](https://garrisondental.com) to see a list of pre-recorded and upcoming live webinars – many with available CE.

Register your email (big red box in the bottom right-hand corner) to get notifications of new CE events delivered right to your inbox!

# Finishing and Polishing Cases

FitStrip™ case



*Images courtesy of Dr. Jeffery Rosenberg, Philadelphia, PA*

eZr™ case



Rally™ case



## TESTIMONIAL

"Interproximal finishing and polishing is a 'must' whenever restoring a Class II, III or IV. However, it is often difficult to achieve. Garrison's new FitStrip™ device is ideal for interproximal trimming (even light orthodontic IPR) and composite finishing."

*Dr. Howard S. Glazer, DDS FAGD  
Fort Lee, NJ*

"I use the eZr™ system on all zirconia and e.max crowns. This product has been fantastic at polishing zirconia...much better than many other brands!"

*Dr. Phillip Howard, DDS  
Seattle, WA*



**F**lexible  
**I**nterproximal finishers with a  
**T**wist


*New*

## Subgingival Trimmers and Finishers!

Subgingival FitStrips are perfectly designed to remove composite overhangs and excess cement quickly and safely from hard-to-reach subgingival areas. These strips are 33% more narrow at only 2.5 mm in height versus standard FitStrips, and with their center, non-abrasive Safety Zone, slip easily between teeth and below the contact area. The medium, fine and super fine durable diamond coating makes quick work of the cleanup and provides a silky-smooth finish.

**Do the Twist!** Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place.

**One handy handle!** FitStrip's simple, attachable handle makes a world of difference for both clinician ease-of-use and patient comfort.



### Stay in the zone

- Non-abrasive Safety Zone protects contact area
- Super-thin – only 0.05 mm!
- Slips easily and safely between teeth

### Narrow focus

- 33% more narrow than standard FitStrips
- Only 2.5 mm in height

**Ready for a marathon!** FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

### Starts bold — finishes clean

1. Medium grit for fast removal
2. Fine grit smooths things out
3. Super fine strips put the perfect finishing touch on things

# FitStrip™

**F**lexible  
**I**nterproximal finishers with a  
**T**wist

## Simple Subgingival Technique:

1. Slide between teeth. Starting with the medium grit (blue), twist the colored barrel until the strip is locked straight. Position the Safety Zone in the affected interproximal space. With gentle apical pressure and a slight back-and-forth motion, slide the strip between the teeth\*. The barrel may be twisted into a curved position at this time to allow the strip to follow the natural tooth contours.
2. Remove the overhang/excess cement. Move Fitstrip™ in a gentle buccal-lingual back-and-forth motion. Remove FitStrip™ and test the area with floss to ensure the desired result. Repeat as necessary.
3. Finish and polish. Place the fine grit FitStrip™ (red) into the interproximal space using the method in step one. Move it in a gentle buccal-lingual back-and-forth motion three or four times. Remove the fine strip and repeat the process with the super fine (yellow) for an ideal polish.

\*If the teeth are bonded together or the contact is too tight to get the 0.05 mm Safety Zone between the teeth, use the FitStrip™ serrated saw (gray barrel) first to gain access.

**NEW!**

## FitStrip™ Subgingival Finishing and Single Sided Kit

FPSK08 contains:

2 handles

10 Strips:

2 subgingival super fine

2 subgingival fine

1 subgingival medium

1 serrated saw

1 each single sided

.08 mm, .10 mm, .13 mm, .18 mm

FitStrip™

**Kits** All kits come in a package of 10 FitStrips with 2 handles.

### NEW Subgingival Finishing Kit

3 subgingival super fine, 3 subgingival fine, 3 subgingival medium,  
1 serrated saw

FPSK07

### NEW Subgingival Finishing and Single Sided Kit

2 subgingival super fine, 2 subgingival fine, 1 subgingival medium,  
1 serrated saw, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

FPSK08



FPSK08

# FitStrip™

**F**lexible  
**I**nterproximal finishers with a  
**T**wist



FPSK05

## IPR Kit

FPSK05 contains:

2 handles

1 IPR gauge

10 Double sided strips (2 of each):

11 mm, .21 mm, .30 mm, 40 mm, .50 mm

## Your Complete Solution for Interproximal Reduction (IPR)!

.40mm and .50mm strips easily accommodate maximum reduction cases

IPR gauges for exact reduction

**IMPROVED!** Attachable handle for greater ease of use and patient comfort

## FitStrip™

**Kits** All kits come in a package of 10 FitStrips with 2 handles.

**Universal Kit** for IPR and finishing/polishing

1 IPR gauge

4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK04

**IPR Kit**

1 IPR gauge

10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK05

**Starter Kit** for finishing/polishing

2 serrated

4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)

4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)

FPSK01

**Single Sided Kit**

2 serrated

8 single sided (2 of each .08 mm, .10 mm, .13 mm, .18 mm)

FPSK02

**Double Sided Kit**

2 serrated

8 double sided (2 of each .11 mm, .15 mm, .21 mm, .30 mm)











FPSK03



**FitStrips** — a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**



#### Refills

		Single Sided	Double Sided	10 pk w/2 handles	4 ct packs (add – 4 to model #)
Super fine 15 Micron		FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick		
Fine 30 Micron		FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick		
Medium 46 Micron		FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick		
Coarse 90 Micron		FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick		
Serrated		FPSSAW 0.05 mm thick			
Medium 46 Micron			FPSTMDS 0.40 mm thick		
Medium 46 Micron			FPSXTMDS 0.50 mm thick		
Super fine 15 Micron		FPSGXFSS 0.08 mm thick	N/A		
Fine 30 Micron		FPSGFSS 0.10 mm thick	N/A		
Medium 46 Micron		FPSGMSS 0.13 mm thick	N/A		

**FitStrip replacement handles** (2-pack) FPSHANDLE-2

IPR Gauge FPSGAUGE



## Mini Polishers for Composite

### Rally gets you to the finish line faster!

The composite finish line that is. One lightning fast step will put a beautiful satin finish on posterior restorations or add the gray high shine polisher for a lustrous gloss on facials and anteriors.

And if you want to be quick off the starting line, use the blue coarse grit for rapid reduction and shaping. Rally Mini Polishers get you across the finish line and on to your next patient in record time!



#### Get to the point – Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.



#### A clean get-a-way!


You don't want to leave a trail of crumbs or spattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.



#### Mini and mighty!


Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.






**Remove –  
Choose the Blue  
(silicon carbide)**

Polishers for fast  
removal and  
contouring



**Smooth –  
Choose the Red  
(aluminum oxide)**

Polishers for  
smoothing and silky  
shine



**Shine –  
Choose the  
Gray (diamond  
impregnated)**

Polishers for ultra  
high shine



*Slow speed right  
angle latch shank.  
Use with water.*

## Rally™

Item	Description	Qty.
FPCCK01	Composite Polisher Complete Kit	60 (10 each)


Item	Refills	Qty.
FPCCK010	Blue coarse grit point composite polisher	30
FPCCK040	Blue coarse grit cup composite polisher	30
FPCCK010	Red medium grit point composite polisher	30
FPCCK040	Red medium grit cup composite polisher	30
FPCCK010	Gray fine grit point composite polisher	15
FPCCK040	Gray fine grit cup composite polisher	15



Intra-oral Adjustment/Polishing  
for Zirconia and Lithium Disilicate

## Adjust and Polish Zirconia and Lithium Disilicate the eZr Way!

### eZr Intra-oral Adjustment:



- **No chatter, chipping or microfracture** – Unlike diamonds bonded directly to a hard metal shaft, eZr's rubber substrate cushions the diamonds to prevent damage to the restoration.
- **No heat build-up** – The carefully engineered substrate does not transfer heat during use.
- **FAST!** – eZr's coarse grit blue grinders make your adjustments fast and convenient.

### eZr Intra-oral Polishing:



- **Unparalleled smooth polish** – Don't take our word for it! – "The eZr system is the only one I've tried that polishes zirconia to smooth!" Todd Brower, DDS, Holland, Michigan.
- **Increased physical stability** – Highly polished zirconia and lithium disilicate is stable and non-abrasive to opposing dentition.
- **No paste required** – Forget the mess and expense of polishing pastes. eZr™ is all you need.
- **No water required** – No need to use water when polishing. Water is only required for high speed.

**eZr™** is a three-step diamond polishing system for the processing of zirconium oxide and lithium disilicate. The grinders and polishers are ideal for adjusting and correcting chairside restorations such as inlays, onlays, veneers and crowns.



*Slow speed right  
angle latch shank  
Steam autoclaveable*



**1**

Start by gross adjusting and contouring using the blue coarse grit grinders. Use the flat end taper for occlusal/labial/lingual/buccal surfaces and the disc for interproximal surfaces.













**2**

Use the green medium grit polishers for smoothing and pre-polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



**3**

Use the orange fine grit polishers for high gloss polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.

Item	Description	Qty.
FPZK01	Complete system kit Gross removal, smoothing, pre-polishing, and high gloss polishing	2 grinders 6 polishers
FPZC030	 Blue coarse grit disc Gross removal	5/pk
FPZC010	 Blue coarse grit flat end taper Gross removal	5/pk
FPZM010 (not in kit)	 Green medium grit small point polisher Smoothing and pre-polishing	5/pk
FPZM020	 Green medium grit large point polisher Smoothing and pre-polishing	5/pk
FPZM030	 Green medium grit disc polisher Smoothing and pre-polishing	5/pk
FPZM040	 Green medium grit cup polisher Smoothing and pre-polishing	5/pk
FPZF010 (not in kit)	 Orange fine grit small point polisher High gloss polishing	5/pk
FPZF020	 Orange fine grit large point polisher High gloss polishing	5/pk
FPZF030	 Orange fine grit disc polisher High gloss polishing	5/pk
FPZF040	 Orange fine grit cup polisher High gloss polishing	5/pk

## Cases using other matrix systems



Images courtesy of Dr. Julio Maya



Images courtesy of Dr. Dory Stutman



Images courtesy of Dr. Matthew A. Nejad, DDS, Beverly Hills, CA



Images courtesy of Dr. Jack D. Griffin, Jr., DMD, MAGD, Weldon Spring, MO



Images courtesy of Dr. Jack D. Griffin, Jr., DMD, MAGD, Weldon Spring, MO

### Reel Matrix™

## TESTIMONIALS

"Using a matrix system like Reel Matrix can enable them (doctors) to quickly and precisely place anatomically correct Class II posterior composite restorations that demonstrate exceptional contacts."

*Dr. Dory Stutman*

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Dr. Matthew A. Nejad, DDS  
Beverly Hills, CA*

"When there is less than 2mm of separation between the adjacent preparations, back-to-back restoring should be considered. In this case, Pinch™ matrices were placed on the first bicuspid and molar and tightened. These feature a blue mylar and an attached tension retainer that is pinched, securing it to the tooth. Slick Band™ sectional matrices are then placed on adjacent restorations and held in place with wedges."


*Dr. Jack D. Griffin, Jr., DMD, MAGD  
Weldon Spring, MO*

"To ensure a tight contact, a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument. This ensures that the contact is tight and well-formed while light curing though it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD  
Weldon Spring, MO*

### PerForm™

## Cases using anteriors/cervicals



Images courtesy of Dr. Zunzarren, DDS



Images courtesy of Dr. Abdi Sameni, DDS


### VariStrip™

How can you produce a natural curvature for anterior composites and eliminate black triangles? Check out Garrison's curved "Varistrips" to recreate the natural interproximal anatomy for anterior teeth.

- Helps avoid incorrect flat embrasures, black triangles and food traps that can be caused by old-fashioned straight mylars.
- No need for wedges on subgingival fillings, as the contours naturally close off the apical aspects for you.
- Excellent for class IVs, diastema closures, etc.
- Variable height allows you to slide each strip in until it precisely matches the tooth height before cutting off excess.
- The blue tint does not affect cure, but ends the annoyance of losing clear strips on patient bibs or bracket tables.

Blue View VariStrips—a must-have super-simple product for every dental office.

*Dr. Michael Curtis, Bridgeport, CT, USA*



### Clinical Tip

After using your preferred bonding procedure, place the restorative material into the preparation in advance of using the cervical matrix.

Place the cervical edge of the matrix against the tooth (Fig. 1).

While pressing the gingival edge of the matrix firmly against the tooth, roll the matrix onto the tooth (Fig. 2). Excess material will be expressed from the occlusal edge of the matrix rather than subgingivally, greatly simplifying clean up.

Cure directly through the matrix.

# ReelMatrix

Real Easy. Real Tight. Real Fast!

## For all posterior restorations!

**1 Real Easy:**  
The reels pop into the handle without any fuss! Press the reel into the handle and retract the grippers to lock the reel in place. Plus, if you use the regular matrices (yellow hubs) they can be easily pushed through contacts.

**2 Real Tight:**  
Tight contacts and excellent contour are what Reel Matrix is all about. These thin matrices are pre-contoured in three dimensions and when combined with anatomical wedges and separator rings, the result is amazingly tight contacts at the natural height of contour.

**3 Real Fast:**  
Pre-loaded reels mean you can just lock them in and go. The retainerless design greatly speeds up quadrant work while improving visibility and patient comfort. Add an FX500 separator ring and you virtually eliminate buccal and lingual flash.



## Margin Elevation Technique restores deep caries without the trauma.

Margin elevation requires a matrix that will properly adapt and isolate the subgingival area. This is not possible with traditional matrix materials (Figure A). The narrow Margin Elevation Band, with its enhanced conical design, is able to slide down subgingivally and properly seal this area (Figure B).




Figure A




Figure B

For the complete margin elevation technique, go to [garrissdental.com](http://garrissdental.com)



**RM500H**  
Margin elevation



**RM100D**  
5mm deadsoft



**RM200D**  
6.3mm deadsoft



**RM100H**  
5mm regular



**RM200H**  
6.3mm regular



**RM100B**  
5mm blue



**RM200B**  
6.3mm blue

The Reel Matrix™ is the ideal delivery system for the Margin Elevation Band. Its extremely small size and light weight prevents the accidental dislodgement that can occur with long, cumbersome metal retainers.



Margin elevation band is .002" thin.  
Super-thin dead-soft bands are .0012" thin.  
Regular bands are .0015" thin and can be placed through existing contacts. A Blue View cure-through option at .0029" thin is available if you prefer trans-enamel curing.

# ReelMatrix

Real Easy. Real Tight. Real Fast!

Item	Description
RMK05	Kit Contains 80 matrices (20 RM500H, 10 ea. RM100D, RM200D, RM100H, RM200H, RM100B, RM200B) 1 tensioning instrument 1 PerForm™ instruments (set of 2)

Refills		Qty.
RM500H	Margin elevation band	50
RM100D	5mm dead-soft matrices	50
RM200D	6.3mm dead-soft matrices	50
RM100H	5mm regular matrices	50
RM200H	6.3mm regular matrices	50
RM100B	5mm blue matrices	50
RM200B	6.3mm blue matrices	50
RMR100	5.9mm empty reel	100
RMR200	6.7mm empty reel	100
RMH-100	Tensioning instrument	

## IMPORTANT ACCESSORIES

PerForm™ Proximal Contact Instrument  
(included in RMK05 kit)

Page 35

Great for core build-ups

# SlickBands™

## Tofflemire-style

**Non-stick technology developed by Garrison®:** Modern bonding agents can adhere to uncoated matrix bands—Slick Bands™ slide right out.

Reaches over to the adjacent tooth more easily to help you create broader, deeper contacts.

- Extremely wide preps are simplified.
- Improved visibility for core buildups.



**Proprietary  
Right Curve™  
design:**

The increased curvature of the all-new **Right-Curve matrix band** produces outstanding gingival seal and reaches over to adjacent teeth easier for broader, deeper contacts. This creates a significantly increased funnel-shaped cone when wrapped around a tooth. It's the right choice for many restorations.


Tighter gingival seal is beneficial for all restorations.

Traditional band with narrow 'funnel'.



## Margin Elevation Band

The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.



A margin elevation matrix will give you the best gingival adaptation and contour for the deep margin areas. It is trimmed in height to allow improved adaptation to deep margin areas.

"We produced the margin elevation band to facilitate the deep margin elevation technique which is becoming increasingly popular with adhesive dentistry. Properly performed, this technique allows for a conservative indirect restoration to be completed without additional invasive treatments. The use of this band and unique matrix shape allows the band to form a tight seal in deeper cervical areas where other bands would usually be unable to reach."

*Matthew A. Nejad, DDS  
Beverly Hills, CA*

# Slick Bands™ Tofflemire-style Matrices


**Clinical Tip:** Use dead-soft green bands when contact is COMPLETELY broken with adjacent teeth for maximum burnishability. When contact is NOT broken with adjacent teeth, use regular gray bands.

## IMPROVED Dead-soft Bands

with non-stick coating



#1 Universal for general restorations




**Right-Curve™** matrices provide a lower contact point and improved isolation. Excellent for all restorations.



#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations




Tofflemire-style matrices are .0014" thick  
Margin elevation matrices are .002" thick

## Regular Bands


with non-stick coating




#1 Universal for general restorations



Right-Curve™ matrices provide a lower contact point and improved isolation. Excellent for all restorations.



#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations



SMTK06

## Slick Bands™

**Description**  
Slick Bands Kit  
300 Assorted Slick Bands  
100 G-Wedges  
PerForm instruments (1 set; shown on page 35)

**Item**  
SMTK06

Dead Soft Bands	Item
#1 Universal Tofflemire	SMT200D10
Right-Curve	SMT400D10
#13 Pedodontic	SMT050D10
#2 Sub-gingival	SMT300D10

Standard (-M)	Bulk
Qty.	Qty.
50	100
50	100
50	100
50	100

Regular Bands	Item
#1 Universal Tofflemire	SMT200H10
Right-Curve	SMT400H10
#13 Pedodontic	SMT050H10
#2 Sub-gingival	SMT300H10

Standard (-M)	Bulk
Qty.	Qty.
50	100
50	100
50	100
50	100

Regular Bands	Item
Margin elevation band	SMT500H10

Standard (-M)	Bulk
Qty.	Qty.
50	100

To order standard packs, add a -M after the item #.

# Pinch<sup>™</sup> Matrices

Your fastest and easiest choice when you need a Tofflemire-style matrix.

Not every Class II can be restored with Composi-Tight®, so when you need a full matrix, Garrison® has you covered. Pinch<sup>™</sup> Matrices are retainerless Tofflemire-style matrices, in cure through or metal, with a built-in tensioning ring. There's no more fiddling with retainers or fancy tightening devices, simply slip around the tooth to be restored and Pinch!

## Retainerless:



The easiest-to-place circumferential matrix eliminates bulky retainers.

## Pre-contoured:

Anatomy-hugging shape helps produce properly contoured restorations.

## Blue Tint:


Provides contrast between matrix and tooth structure without compromising composite polymerization.



## Pinch Matrices


Metal Kit	Description	Qty.
KMMK1	Metal Assortment Kit 50/ea. size	150
<b>Metal Refills</b>		
KMM100	Metal Premolar	50
KMM175	Metal Molar	50
KMM200	Metal Tall Molar	50
<b>Blue View Kit</b>		
KPMK1	Cure Through Assortment Kit 50/ea. size	100
<b>Blue View Refills</b>		
KPM100	Cure Through Premolar	50
KPM200	Cure Through Molar	50

## How PerForm™ works:




1A

**STEP 1:**  
Place the first composite increment. Press PerForm™ into the uncured composite increment and either push (distal; fig. 1A) or pull (mesial; fig. 1B).



1B

**STEP 2:** Place the curing light directly on the Fresnel lens and cure.



2

# PerForm

Proximal Contact Instruments

## PerForm puts you in control.

- Improves contacts with any matrix system
- Control the height and breadth of contacts
- Increase the tightness of contacts
- Unique lens controls curing light diffusion
- Steam autoclave only

### PerForm™

Item	Kit	Qty.
CFK01	PerForm (1/ea. size)	2/pk
CF100	1.1 mm Small Instruments	2/pk
CF200	2 mm Large Instruments	2/pk



"To ensure a tight contact a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument. This ensures that the contact is tight and well-formed while light curing through it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD  
Weldon Spring, MO*




Photo  
courtesy of  
David S. Hornbrook, DDS, FAACD  
International Lecturer and  
Private Practice, San Diego, CA

**Curved Form/tapered:** Provides the optimum curvature and band height for virtually every anterior restoration.

**Blue Tint:** Provides contrast between matrix and tooth structure without compromising composite polymerization.

## Our Most Popular Anterior Matrix

**BlueView™**  
**VariStrip™**  
Contoured Anterior Matrices



The Blue View™ VariStrip™ is a revolutionary contoured anterior matrix providing the optimum curvature and band height for virtually every anterior restoration. The 0.0020" thin plastic anatomical strip is tapered from one end to the other. You can place the strip interproximally and then slide it until the tooth height is perfectly matched. The pre-contouring easily recreates the occluso-gingival anatomy, avoiding flat embrasures.

Blue View™ VariStrip™

Item	Kit	Qty.
VS01	VariStrip	100

# BlueView™

## Cervical Matrices

"I love the Blue View Cervical Matrix! It works so well. My assistants love watching me use this system. The patients love the excitement in our faces. It makes the process easy and fast. Gotta love that! We are also hooked on using the ReelMatrix. Thanks so much for making great products."

Cynthia Brattesani, DDS,  
San Francisco, CA, USA



### Flexible

- adapts to individual anatomy for improved contours
- use it to directly compact the composite into the restoration
- thin edges allow excess composite to flow out for improved anatomy

### Blue Tint


- provides contrast between matrix and tooth structure
- does not compromise composite polymerization

### Convenient!

- protects the restoration from contamination
- eliminates time spent hand sculpting
- eliminates air-inhibited layer
- reduces finishing time
- 7 sizes for a perfect fit
- Gingival Retractor included!

### KCMK1 INCLUDES

- 275 Assorted Matrices
- 1 Matrix placement Instrument
- 1 Gingival Retractor



### MATRICES INCLUDED:

- Anterior/Premolar  
KCMA60 - 50 pcs
- Molar  
KCMA70 - 50 pcs
- KCMA10 - 35 pcs
- KCMA20 - 35 pcs
- KCMA30 - 35 pcs
- KCMA40 - 35 pcs
- KCMA50 - 35 pcs

## Blue View™ Cervical Matrices

### Item

KCMK1

### Description

Assortment Kit w/Instruments

### Refills

KCMA10

— Refill 100 pack

KCMA20

Refill 100 pack —

KCMA30

— Refill 100 pack

KCMA40

Refill 100 pack —

KCMA50

— Refill 100 pack

KCMA60


Refill 100 pack —

KCMA70

— Refill 100 pack

KCMINST1

Placement Instrument




# Composi-Tight® Matrix Forceps

Precise Placement, Unbeatable Grip.



BFA-1



**Precise Placement:** The tips have a 75-degree angle with an occlusal-gingival orientation, providing a more appropriate placement angle than cotton forceps. The petite, carbide tips allow the matrix to be grasped at its upper edge to prevent distortion of the interproximal area of the band.

**Unbeatable Grip:** The tungsten carbide tips provide exceptional grip and tactile control for tough access cases.



## Check out the GRIP!

Effortlessly remove the matrix bands that otherwise would be a struggle.

### Matrix Forceps

BFA-1

Angled tip for occlusal placement

# Composi-Tight® 3D Fusion™

Composi-Tight® 3D Fusion™ includes dentistry's premier heavy-duty stainless steel ring placement forceps. Regardless of your hand size or strength, Composi-Tight® 3D Fusion™ forceps will give you incredible control for the precise placement of super-strong separator rings over even the widest molars.


**Composi-Tight® 3D Fusion™  
Ring Placement Forceps**  
FXP01

Composi-Tight® 3D Fusion™ forceps will work with virtually any ring from any manufacturer.



## Ring Placement Forceps

Holds all your sectional matrix system rings.



**Ring Placement Forceps**  
AUMRDF-100

**Universal Forceps**  
for placing rubber dam  
clamps and rings  
UNVF-100



# Composite Instruments Nonstick titanium-nitride coated

Our multi-function instruments have double condenser configurations and carving blades on opposite ends. Now you can place, pack, shape and carve your way to perfectly contoured composite restorations without ever having to switch instruments!

- Large-diameter, lightweight handles for improved ergonomics.
- Provides non-stick surface for placement and shaping of composites.
- 40% smoother than stainless steel for a non-porous, non-stick surface.
- Withstands diverse cleaning/sterilization methods, including ultrasonic cleaning.



“Skippy Special”  
TN008

Item	Description
TN008	“Skippy Special” 2 condensers, ball burnisher, blade, acorn



5-in-1  
TN009

Item	Description
TN009	5-in-1 Multi-Function Instrument 2 blades, 2 condensers, acorn

**Item Description**

**TN001 Medium Blades**  
 Placement of composite material  
 Gingival retraction to sculpt under free margins  
 Burnishing matrix bands



**TN002 Medium Blade with Medium Plugger**  
 Contouring small Class I, II, V  
 Packing composite material into conservative preparations  
 Burnishing matrix bands



**TN003 Medium Blade with Large Plugger**  
 Contouring large Class I, II, V  
 Packing composite material into large preparations  
 Convex plugger tip minimizes pullback of composite material



**Item Description**

**TN004 Extra Thin Blades**  
 Reach under margins to hone and sculpt  
 Sculpting interproximally  
 Packing cord in gingival crevice



**TN005 Long Extra Thin Blades**  
 Allows additional access for hard-to-reach  
 posterior contouring and shaping  
 Sculpting interproximally  
 Packing cord in gingival crevice








**TN006 Small/Medium Multi-Function Plugger**  
 Placement of composite material in smaller Class II and III  
 Packing material toward gingival seats  
 Carving occlusal anatomy with acorn carver



# Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials	Page	Special Notes
<b>Distal of canine, short or malpositioned</b> 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! All new separator rings designed to hold firmly on tough cases.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
<b>Conservative Class II</b> 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
<b>Slightly wide and/or deep</b> 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.

Procedure	Recommended Materials	Page	Special Notes
<b>Large Class II with missing cusps</b> 	Composi-Tight® 3D Fusion™ Sectional Matrix Kit—Green Wide Prep Ring	4–5	EXCLUSIVE! Dentistry's first Wide Prep ring greatly simplifies these more challenging restorations.
	Reel Matrix™—Metal or Blue View	30–31	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Pinch™ Matrices—Metal or Blue View	34	Very fast retainerless Tofflemire-style matrices.
	Slick Bands™ Tofflemire-Style Matrices	32–33	Traditional Tofflemire-style with non-stick coating.
<b>Small to medium size Class III</b> 	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
<b>Large Class III or Class IV</b> 	VariStrip™ Contoured Anterior Matrices	36	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
<b>Class V</b> 	Blue View™ Cervical Matrix	37	Contours and directly compresses resin while eliminating air-inhibited layer.
<b>Crown build-ups, post and core</b> 	Reel Matrix™	30–31	Best choice. Margin Elevation band available for very deep buildups. Retainerless convenience.
	Pinch™ Matrices—Metal or Blue View	34	Very fast; may not provide as tight a cervical seal as other matrices.
	Slick Bands™ Tofflemire-Style Non-Stick Matrices	32–33	Margin Elevation band available for very deep buildups if you prefer a traditional matrix retainer.

# Order Form

Fax: 616.842.2430

Item Number	Description	Quantity	Price/Unit	Total Price
		Shipping		
		Sales Tax		
		Total Price		

Call for additional express shipping rates and delivery. All prices subject to change without notice. Appropriate shipping charges and sales tax will be added. Special offers not valid with any other advertised specials.

Doctor's Name		
Contact Name		
Address		
City	State	Zip
Phone	Fax	
Email		
Payment Method		
Credit Card: <input type="checkbox"/> Visa <input type="checkbox"/> MasterCard <input type="checkbox"/> Am. Express <input type="checkbox"/> Discover <input type="checkbox"/> C.O.D. (add \$13.00)		
Credit Card #	CC Ref #	
Exp. Date	Signature	
Credit Card Billing Information		
Address		
City	State	Zip

toll free 888.437.0032 | [garrisondental.com](http://garrisondental.com)

# Customer Service Information

Call **888.437.0032 Toll Free**  
(USA & Canada)  
616.842.2244

Fax 616.842.2430

Email [gds@garrisondental.com](mailto:gds@garrisondental.com)

Online **[garrisondental.com](http://garrisondental.com)**  
Look for our Learning Center with clinical tips, case studies and technique guides.

Mail 150 DeWitt Lane  
Spring Lake, MI 49456  
USA

# Shipping Rates

Regular UPS	2nd Day Air	Next Day Air	C.O.D.
\$9.75	\$19.00	\$34.00	+\$13.00

Express orders must be placed by 4pm EST to ship same day.

Shipping charges subject to change without notice.

# Payment Methods

VISA, MasterCard, American Express, Discover, Check Prepay or C.O.D.

Prices are subject to change without notice.

# Distribution

Garrison is a global leader for practical and innovative dental products and accessories, with distributors worldwide.

We are pleased to announce that your favorite Garrison products are available through several full-service dealers in the USA.

To find a dealer, contact our headquarters in Spring Lake, Michigan. Phone: 616.842.2244 or Email: [gds@garrisondental.com](mailto:gds@garrisondental.com)



Order Toll-Free • Order by Fax • Order Online • Order by Mail

## Our Guarantee – It's Simple!

### **Satisfaction = 100%**

Everything you order from Garrison® Dental Solutions is backed by our risk-free six-month money-back **Satisfaction Guarantee**. If you are not satisfied with your product, return the unused portion and we will issue a refund.

### **Quality = 100%**

All Garrison products are covered by our **Quality Guarantee**.

We will repair or replace any Garrison® product that is shown to be defective in materials or workmanship based on the individual product's anticipated lifetime. Please visit our website for details on our various products, [garrisondental.com](http://garrisondental.com).