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Garrison®

English

This cleaning, disinfection and sterilization guide has been provided in accordance with ISO 17664, covering the proper techniques necessary to ensure Garrison® Dental Solutions' product will be suitable for safe and optimal use.

Failure to follow the stated procedure may result in product failure, premature wear and/or non-optimal performance. Any use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.

Contact Information:

For the most current information, visit www.garrisondental.com/sterilization

To contact the manufacturer directly:

U.S./Canada:
Garrison® Dental Solutions, LLC
150 DeWitt Lane
Spring Lake, MI 49456 USA
Tel: +1.616.842.2244
Fax: +1.616.842.2430
Email: gds@garrisondental.com

Europe Office:
Garrison® Dental Solutions
Carlstrasse 50
D-52531 Uebach-Palenberg
Germany
Tel: +49.(0)2451.971.409
Fax: +49.(0)2451.971.410

How to use this guide:

Identify product requiring cleaning, disinfection or sterilization in the left column of the product matrix. Follow the procedure indicated in the right column of the product matrix.

Product Line	Cleaning, Disinfection & Sterilization Procedure:
Sectional Matrix Kits	A, B, C, E, F (as appropriate to kit contents)
Matrix Ring: SG400, SG500, SG600, AU400, AU500, G100, G200, FX400, FX500, FX600, 3DXR, 3D500, NT400, NT500, NTG400, NTG500, Q400L-1, Q400R-1, Q500L-1, Q500R-1	A
Matrix Forceps – All varieties	A
Ring Forceps – All varieties	A
TN Series Composite Instruments	A
FitStrip™ Interproximal Finishers	A
PerForm™ Instruments	A
Blue View™ Cervical Handle	A
Matrix Bands: SG Series, AN Series, FX Series, FXH Series, SXR Series, M Series, SM Series, AU Series, B Series, NTG Series, NTM Series, SMT Series	B
Wedges: Strata-G™, Fusion™ Anterior, Composi-Tight® 3D Fusion™, Wedge Wands®, G-Wedge™, FenderWedge™, NiTin Wedges ,Quad™ Wedge (includes assortment kits)	B
Blue View™ VariStrip™	B
ReelMatrix™ Spools & Matrices	C
ReelMatrix™ Instrument	C
eZr™ Polishing/Adjusting Instruments	D
A+Wedge® (includes assortment kits)	E
Rally™ Polishers	F
ProxiBur™ Interproximal Reduction Burs	G

*not a trademark of Garrison® Dental Solutions, LLC

Procedure "A"

Basic Principles:

- All products are supplied non-sterile and may be appropriately cleaned, disinfected, and sterilized before use.
- Ensure equipment used is serviced and checked regularly, and manufacturer's use instructions are followed.
- Please heed regulations applicable to your country.
- Repeated reprocessing may affect the product. Decommissioning of the device is normally determined by wear and damage prior to use.
- Avoid contact of dissimilar metals at any point during the cleaning & sterilization process.

For users employing automated cleaning or disinfection equipment:

- Garrison Dental Solutions' reusable devices are intended for terminal steam sterilization as the validated microbial inactivation step. Automated cleaning or disinfection equipment may be used prior to sterilization, provided it does not alter device materials or intended use. However,no specific automated disinfection process is claimed or validated for these devices.
- The use of automated equipment should always be followed by validated manual cleaning and terminal steam sterilization as described in the Instructions for Use. No additional disinfection validation is required or endorsed under current product claims.
- Any use of other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.

Cleaning:

- Products should be cleaned as soon as possible after use to prevent soils from hardening.
- Remove all visible debris after use and prior to sterilization.
- Product can be manually cleaned with IMS® Daily Clean*, following manufacturer's instructions.
- Use a nylon brush or similar soft brush. For SG400, SG500, SG600, FX400, FX500, FX600, 3D500, 3DXR, Q400L-1, Q400R-1, Q500L-1, & Q500R-1: Take care not to damage the Soft-Face™ during cleaning.
- You may optionally pre-soak the contaminated rings with IMS® Daily Clean*, following the manufacturer's instructions, for up to 6 minutes to aid in manual cleaning.
- DO NOT expose product to solutions containing chemicals such as alcohol, acetone, ketone, or other harsh chemical disinfectants.
- Visually inspect product for damage (eg. rust and cracking). Remove from service if damage is present.
- If any visible debris remains, repeat process above as needed.

Sterilization:

- Process product through steam autoclave and store only in CE/FDA approved sterilization pouch.
- Steam autoclave sterilization only. Perform one cycle at a minimum of 270°F (132°C) and 30.46 psi (210kPa) for a minimium of 4 minutes, following the autoclave manufacturer's instructions, with a 20–30 minute drying time. Never exceed temperatures of 278.6°F (137°C).
- DO NOT use chemical sterilization techniques.

Storage:

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

Procedure "B"

Basic Principles:

- All products are supplied non-sterile and may be appropriately sterilized before use.
- Ensure equipment used is serviced and checked regularly, and manufacturer's use instructions are followed.
- Please heed regulations applicable to your country.
- Products are single use only and cannot be reprocessed.
- Avoid contact of dissimilar metals at any point during the sterilization process.

- Any use of other reprocessing procedures/methods is at the discretion and sole responsibility of th epractitioner.

Sterilization:

- Process product through steam autoclave and store only in CE/FDA approved sterilization pouch.
- Steam autoclave sterilization only. Perform one cycle at a minimum of 270°F (132°C) and 30.46 psi (210kPa) for a minimum of 4 minutes, following the autoclave manufacturer's instructions, with a 20–30 minute drying time. Never exceed temperatures of 278.6°F (137°C).
- DO NOT use chemical sterilization techniques.

Storage:

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

Procedure "C"

Basic Principles:

- Ensure equipment used is serviced and checked regularly, and manufacturer's use instructions are followed.
- Please heed regulations applicable to your country.
- Avoid contact of dissimilar metals at any point during the sterilization process.
- Any use of other reprocessing procedures/methods is at the discretion and sole responsibility of th epractitioner.

Spools and Matrices Cleaning and Sterilization:

- Spools and Matrices are single use disposable items. Do not attempt to sterilize them. Prior to use they may be sprayed with a chemical diinfectant if desired.

Instrument Cleaning:

- Instrument must be in open position. Take care not to damage tips.
- Remove visible debris or organic matter from the instrument and surface with a disposable wipe or surface brush. This is to be done by scrubbing with warm tap water and soap or detergent.
- Inspect instrument for debris. Repeat cleaning if there is debris left.
- All external surfaces must be wiped with a lint-free cloth or sprayed with a suitable disinfectant.
- Rinse external surfaces with warm tap water.
- Dry instrument.
- Do not immerse the instrument in liquid or place it in an ultrasonic cleaning unit.

Instrument Sterilization:

- Sterilize in a steam autoclave using distilled water at a temperature of 134°C (273°F) for 3 minutes.
- Sterilization is NOT a substitute for proper cleaning. Clean the instrument first then sterilize it.
- Do not use cold sterilization methods that involve soaking the instrument in a liquid.

Instrument Storage:

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

Procedure "D"

Basic Principles:

- All products are supplied non-sterile and may be appropriately cleaned, disinfected and sterilized before use.
- Ensure equipment used is serviced and checked regularly, and manufacturer's use instructions are followed.
- Please heed regulations applicable to your country.
- Repeated reprocessing may affect the product. Decommissioning of the device is normally determined by inspection for wear and damage prior to use.
- Avoid contact of dissimilar metals at any point during the cleaning & sterilization process to prevent corrosion.

Notices:

- Avoid the use of strong acids and bases as they may oxidize the stainless steel shaft.
- Rinse with distilled water after treatment with cleaning and disinfectant solutions.
- Avoid temperatures >150°C (302°F).
- Ultrasonic bath temperatures must not exceed temperatures of 42°C (107°F) to avoid possible coagulation of proteins resulting from blood and tissue.
- The label on the cleaning and/or disinfecting solution MUST specifically say "suitable for rubber polishers or synthetics/silicones."

Point of Use:

- Products should be cleaned as soon as possible after use to prevent soils from hardening. Pre-clean under running water with a plastic bristle brush.

Manual Cleaning:

- Under running water with a plastic bristle brush.
- Optional – Ultrasonic assisted cleaning with a suitable cleaning agent and disinfectant. Rinse afterwards under running water.

Automatic Cleaning:

- Thermal disinfector – Manufacturer's specification according to EN ISO 15883 using the cleaning program as indicated by the manufacturer's operating instructions.

Disinfection:

- Use a solution classified suitable for rubber and silicone polishers/synthetics by the manufacturer. Follow manufacturer's recommendations for exposure times and concentrations.
- Rinse with distilled water after treatment with cleaning and disinfectant solutions.

Drying:

- Dry with fresh, clean, lint-free cellulose tissues.

Inspection:

- Visually check all instruments for cleanliness and safe suitability of use with 5–10X magnification. Repeat cleaning if necessary. Decommission and dispose of instruments in the event of discernible defects.

Sterilization:

- Steam Autoclave according to manufacturer's operating instructions using a hold time of 5 minutes full cycle, a sterilization temperature of 134°C (273°F) and a drying time of 10 minutes.

Storage:

- Store instruments packed and protected in suitable sterile packaging, cassettes or retainers in dry, clean conditions at ambient temperature.

Validated Materials:

Manual Cleaning – Durr Dental ID 212, 2% concentration, 5 minute holding time
Ultrasonic – Durr Dental ID 212 ready made solution, 2% concentration, 2 minute holding time
Automatic Cleaning – Miele G 7883, Special Program 93°C (199°F) for 10 minutes, 10 minute holding time, detergent – Dr. Weigert neodisher MediClean Dental, rinsing for 3 minutes at 75°C (167°F) with neodisher Z Dental
Sterilization – Holding time full cycle 5 minutes, sterilization temperature 134°C (273°F), drying time 10 minutes, threshold values of contents for feed-water and steam condensates

Procedure "E"

Basic Principles:

- Products are single use only and cannot be reprocessed.

Sterilization:

- DO NOT clean or sterilize as this will affect the astringent coating.

Storage:

- Store in original packaging in a cool dry location.
- Do not mix in with similar products.

Procedure "F"

Basic Principles:

- All products are supplied non-sterile and may be appropriately sterilized before use.
- Ensure equipment used is serviced and checked regularly, and manufacturer's use instructions are followed.
- Please heed regulations applicable to your country.
- Products are single use only.

Sterilization:

- Moist heat sterilization with fractional prevacuum. Full cycle, 134°C (273°F) / 3.1 bar (45 psi) for at least 5 minutes, drying time 10 minutes. Do not exceed 135°C (275°F). IMPORTANT: Use suitable packaging for sterilization.
- Warning: not suitable for sterilization with hot air or chemiclave. This procedure was validated by an independent and accredited laboratory.

Storage:

- A dry, dust-free storage at 2–28°C (36–82°F) is recommended. Improper storage conditions can shorten the shelf life and cause malfunction of the product.

Procedure "G"

Basic Principles:
Cleaning and sterilization parameters referenced below are based on validation performed by Strauss & Co., the contract manufacturer.

- Garrison's dental Diamond burs are designed for several uses and are eligible for both cleaning and autoclave sterilization.
- Clean and sterilize the burs before first use and before each reuse.
- Use only a neutral pH cleaning / disinfecting solution which is approved for this use (through CE Marking or FDA approval/clearance) and compatible with Stainless-Steel and Nickel.
- The cleaning / disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor.
- Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments.
- Cleaning may be performed using manual cleaning or automated cleaning. Sterilization should be conducted using steam sterilization, per the instructions below for each method.

Manual Cleaning:

- Immediately after use, immerse all instruments in a disinfection solution compatible with stainless steel and Nickel, according to the cleaning agent manufacturer's instructions (this product has been validated by the Manufacturer using Cidex OPA in concentration of 0.3%). Use the lowest recommended concentration and lowest recommended temperature (an excessive concentration may cause corrosion or other defects on instruments).
- Immerse the products for 5 minutes.
- While soaked scrub the device thoroughly with a soft nylon brush for 5 minutes.
- Allow drying for at least 10 minutes. Make sure that the product is completely dry.
- Observe the device to assure that no discoloration or corrosion signs or other visual contamination appeared.
- Move forward to the steam sterilization process.

Automated Cleaning:

- As ***manual pre-cleaning***, soak all instruments immediately after use in a detergent and disinfecting neutral pH solution combined with proteolytic enzyme for 5 minutes or as specified by the cleaning agent manufacturer. Use the lowest recommended concentration by the manufacture. Endozime® Premium with APA (Advanced Proteolytic Action) detergent was utilized for this product's validation by the Manufacturer at 0.1% equal to 1ml/liter in lukewarm water.
- Follow instructions and observe concentrations and immersion times given by the specific cleaning agent's manufacturer (an excessive concentration may cause corrosion or other defects on instruments).
- Under running water, scrub the device thoroughly with a soft nylon brush to remove any remaining blood or debris.
- Rinse under tap water for one (1) minute.
- Place the burs in the automatic washer (until they are placed in the washer avoid contact between them.
- **Utilize the automated process:** Run the automatic wash cycle according to automated washer manufacturer's instructions.
- **The following parameters were utilized by the contract manufacturer (Strauss & Co.) to validate this product:**
 - Two (2) Minutes prewash at 30±5°C (86±41°F) with tap water.
 - Ten (10) Minutes main wash with the neutral pH detergent and tap water at 113±41°F (45±5°C).
 - One (1) Minute rinse with tap water at 86°F (30°C).
 - Ten (10) Minutes distilled water rinse at 86°F (30°C).
 - Twenty (20) Minutes air drying phase at high temperature of 212°F (100°C).
- Observe the device to assure that no discoloration or corrosion signs or other visual contamination appeared.
- Move forward to the steam sterilization process.

Steam Sterilization Instructions:

- Follow the autoclave manufacturer's instruction to sterilize the products. Do not exceed the maximal load allowed by the manufacturer.
- Note: Local infection control practice may recommend a different combination of holding time and temperature.
- Place the products into the steam sterilizer into a sterile barrier pouch.
- Note: Use microbial barrier pouch that is cleared by the US FDA (US users) and CE-Marked approved (EU users).
- Use an autoclave that meets applicable regulatory requirements for the intended market (e.g., FDA-cleared for U.S., CE-marked for EU).
- Note: Consult the autoclave's manufacturer instructions for use.
- Method validated by the Manufacturer:
 - Wrapped instruments, Gravity Displacement.
 - Temperature: 275^0F/ 135^0C.
 - Exposure Time: 10 minutes.
 - Drying Time: 30 minutes.
- After drying, make sure the product is completely dry.
- Use maximal loaded autoclave with wrapped products.
- After sterilizing the device, do not remove burs from sterile pouch until used. If burs are being sterilized unwrapped, use them immediately after sterilization. Use aseptic technique.

Storage:

- After sterilization, keep the devices in sterilization packaging in a dry and clean environment, under room temperature.
- For additional information, refer to the Strauss & Co. Instructions for Use for reusable dental burs.

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