


Garrison BAND COMPARISON CHART

Dental Solutions



FX Series
SXR Series
M – SM – AU Series
B Series

	Pediatric	Bicuspsids	Bicuspsids with deep cervical preps	Bicuspsids and small molars	Molars	Molars with deep cervical preps						
FXH Series	Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.											
	thickness 0.0015"											
	height	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm					
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								FXH100	50	\$78	100	\$137
								FXH150	50	\$78	100	\$137
								FXH175	50	\$78	100	\$137
								FXH200	50	\$78	100	\$137
								FXH300	50	\$78	100	\$137
FX Series	Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.											
	thickness 0.0016"											
	height	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm					
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								FX100	50	\$81	100	\$142
								FX150	30	\$81	60	\$142
								FX175	50	\$81	100	\$142
								FX200	50	\$81	100	\$142
								FX300	30	\$81	60	\$142
NTG Series	NiTin™ full contour matrix bands for ideal anatomy.											
	thickness 0.0015"											
	height	NA	NTG100 4.0 mm w/ext 4.4 mm	NTG150 4.1 mm w/ext 6.0 mm	NTG175 5.0 mm w/ext 5.6 mm	NTG200 5.4 mm w/ext 6.6 mm	NTG300 6.0 mm w/ext 8.7 mm					
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								NTG100	25	\$36	50	\$64
								NTG150	25	\$36	50	\$64
								NTG175	25	\$36	50	\$64
								NTG200	25	\$36	50	\$64
								NTG300	25	\$36	50	\$64
SXR Series	Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.											
	thickness 0.0016"											
	height	NA	SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm w/ext 8.6 mm					
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								SXR100	50	\$73	100	\$127
								SXR150	30	\$73	60	\$127
								SXR175	50	\$73	100	\$127
								SXR200	50	\$73	100	\$127
								SXR300	30	\$73	60	\$127
SM Series	Slick Bands™ Matrices have non-stick coating.											
	thickness 0.0016"											
	height		SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm				
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								SM050			60	\$127
								SM100	50	\$73	100	\$127
								SM150	30	\$73	60	\$127
								SM175	50	\$73	100	\$127
								SM200	50	\$73	100	\$127
								SM300	30	\$73	60	\$127
M Series	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.											
	thickness 0.0015"											
	height		M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm				
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								M050	25	\$60	50	\$92
								M100	50	\$60	100	\$92
								M150	25	\$60	50	\$92
								M175	50	\$60	100	\$92
								M200	50	\$60	100	\$92
								M300	25	\$60	50	\$92
AU Series	The Gold bands are the same as the Composi-Tight® bands.											
	thickness 0.0015"											
	height		AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm				
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								AU050	25	\$60	50	\$92
								AU100	50	\$60	100	\$92
								AU150	25	\$60	50	\$92
								AU200	50	\$60	100	\$92
								AU300	25	\$60	50	\$92
B Series	The original bands have a slightly shorter length.											
	thickness 0.0015"											
	height		B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm				
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								B050	25	\$52	50	\$79
								B100	50	\$52	100	\$79
								B150	25	\$52	50	\$79
								B200	50	\$52	100	\$79
								B300	25	\$52	50	\$79
MC Series	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.											
	thickness 0.0020"											
	height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA					
								Standard (-M)	Bulk			
								Qty.	Price	Qty.		
										Price		
								MC100	50	\$73	100	\$127
								MC175	50	\$73	100	\$127
								MC200	50	\$73	100	\$127

To order standard packs, add a -M after the item #.

© 2021 Garrison Dental Solutions, LLC