

New and Improved Slick Bands™ Tofflemire-style Matrix Kit

When a sectional matrix just won't do: The Slick Bands Tofflemire-style Matrix Kit will!



- **Non-stick technology** developed by Garrison
Modern bonding agents can adhere to uncoated matrix bands — Slick Bands slide right out.
- **NEW! Improved tear-resistance**
Only .0014" thick, the dead-soft green bands are made from a special blend of stainless steel that resists tearing yet is easily burnished.
- **Exclusive Garrison shapes**

Right Curve™

Reaches over to the adjacent tooth more easily to help you create broader, deeper contacts.


- Extremely wide preps are simplified.
- Improved visibility for core buildups.



Tighter gingival seal is beneficial for all restorations.

All new Margin Elevation Band!

The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.



Slick Bands™ Tofflemire-style Matrix Kit

Color-coded
and conveniently organized




The **margin elevation matrix** is trimmed in height to allow proper adaptation and seal of these deep areas so that a composite margin elevation can be performed.

Dead-soft Bands with non-stick coating



#1 Universal for general restorations



Right-Curve™ matrices provide a lower contact point and improved isolation. Excellent for all restorations.




#13 Pedodontic for short and primary teeth



#2 Sub-gingival for deep restorations

Tofflemire-style matrices are .0014" thick
Margin elevation matrices are .002" thick

Regular Bands with non-stick coating



Slick Bands™

Description	Item
Slick Bands Kit	SMTK06
300 Assorted Slick Bands	
100 G-Wedges	
PerForm instruments (1 set)	

Dead Soft Bands	Item	Standard Qty.
#1 Universal Tofflemire	SMT200D10-M	50
Right-Curve	SMT400D10-M	50
#13 Pedodontic	SMT050D10-M	50
#2 Sub-gingival	SMT300D10-M	50

Dead Soft Bands	Item	Bulk Qty.
#1 Universal Tofflemire	SMT200D10	100
Right-Curve	SMT400D10	100
#13 Pedodontic	SMT050D10	100
#2 Sub-gingival	SMT300D10	100

Regular Bands	Item	Standard Qty.
#1 Universal Tofflemire	SMT200H10-M	50
Right-Curve	SMT400H10-M	50
#13 Pedodontic	SMT050H10-M	50
#2 Sub-gingival	SMT300H10-M	50

Regular Bands	Item	Bulk Qty.
#1 Universal Tofflemire	SMT200H10	100
Right-Curve	SMT400H10	100
#13 Pedodontic	SMT050H10	100
#2 Sub-gingival	SMT300H10	100

Regular Bands	Item	Standard Qty.
Margin elevation band	SMT500H10-M	50

Regular Bands	Item	Bulk Qty.
Margin elevation band	SMT500H10	100

PerForm


Description	Item	Qty.
Perform	CFK01	1/ea. size
1.1 mm Small	CF100	2/pk
2 mm Large	CF200	2/pk

Kit includes PerForm™ Proximal Contact Instruments.

They're the fast and easy way to give you the contacts you demand!



- Control the height and breadth of contacts
- Increase the tightness of contacts
- Unique lens controls curing light diffusion



For more information go to www.garrisdental.com. Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

To order: Call toll free **888.437.0032**
or online at www.garrisdental.com



Mention Code
IN0216TF