


Band Comparison Chart

FX Series
SXR Series
M – SM – AU Series
B Series



	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps	Regular Item	Bulk Item
FX Series	3D Fusion™ 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.						FX100-M FX150-M FX175-M FX200-M FX300-M	FX100 FX150 FX175 FX200 FX300
	thickness 0.0016" height	NA	FX100 4.4 mm	FX150 6.0 mm	FX175 5.6 mm	FX200 6.6 mm	FX300 8.7 mm	
SXR Series	Slick Bands™ XR Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.						SXR100-M SXR175-M SXR200-M SXR150-M SXR300-M	SXR100 SXR175 SXR200 SXR150 SXR300
	thickness 0.0016" height	NA	SXR100 4.6 mm	SXR150 3.8 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm	
B Series	Composi-Tight® Original The original bands have a slightly shorter length.						B050-M B100-M B150-M B200-M B300-M	B050 B100 B150 B200 B300
	thickness 0.0015" height		B050 3.2 mm	B100 4.6 mm	B150 3.8 mm	NA	B200 6.4 mm	B300 6.4 mm
SM Series	Slick Bands™ Slick Bands™ Matrices have non-stick coating.						SM100-M SM150-M SM175-M SM200-M SM300-M	SM100 SM150 SM175 SM200 SM300
	thickness 0.0016" height	NA	SM100 4.6 mm	SM150 3.8 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm	
M Series	Composi-Tight® Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.						M050-M M100-M M150-M M175-M M200-M M300-M	M050 M100 M150 M175 M200 M300
	thickness 0.0015" height		M050 3.2 mm	M100 4.6 mm	M150 3.8 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm
AU Series	Composi-Tight® Gold® The Gold bands are the same as the Composi-Tight® bands.						AU050-M AU100-M AU150-M AU200-M AU300-M	AU050 AU100 AU150 AU200 AU300
	thickness 0.0015" height		AU050 3.2 mm	AU100 4.6 mm	AU150 3.8 mm	NA	AU200 6.4 mm	AU300 6.4 mm
MC Series	Composi-Tight® 3D Clear Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.						MC100-M MC175-M MC200-M	MC100 MC175 MC200
	thickness 0.0020" height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA	