


Composi-Tight™ Matrix Forceps

Precise Placement, Unbeatable Grip.



Check out the GRIP! Order a new Matrix Forceps for your office and put them to the test. (Remember, there's a 6-month money-back guarantee!) Try a simple tug-of-war contest with your current instrument. If you are not completely convinced that the Matrix Forceps are a better way to place and remove matrix bands, send them back for a full refund.



Precise Placement: The tips have a 75-degree angle with an occlusal-gingival orientation, providing a more appropriate placement angle than cotton forceps. The petite, carbide tips allow the matrix to be grasped at its upper edge to prevent distortion of the interproximal area of the band.


Unbeatable Grip: The tungsten carbide tips provide exceptional grip and tactile control for tough access cases.

Matrix Forceps
BFA-1 Angled tip for occlusal placement

Ring Placement Forceps

Holds all your sectional matrix system rings.


Updated grip to easily spread rings with less hand fatigue, even over large molars.



Ring Placement Forceps
AUMRDF-100 Ring placement forceps

Universal Forceps

Holds all your rubber dam clamps and sectional matrix retainer rings.



Universal Forceps
UNVF-100 Universal forceps

For more information go to www.garrisondental.com. Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

To order:
Call toll free **888.437.0032**
or online at www.garrisondental.com



Mention Code
IN0216TN