

Rally™

Mini Polishers for Composite

Rally gets you to the finish line faster!

The composite finish line that is. With zero reprocessing and sterilization hassles and a lustrous diamond final gloss, Rally mini-polishers get you across the finish line and on to your next patient.

Get to the point – Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.

Mini, but mighty!

Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.

A clean get-a-way!

You don't want to leave a trail of crumbs or scattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.

Give your mechanic the day off! (Maybe your assistant too!)



Rally is fully assembled – no little mandrels to lose track of, no discs to snap on, there's just nothing you need to mess with – including cleaning and sterilization.



Garrison
Dental Solutions


Rally™

Mini Polishers for Composite




Remove - Choose the Blue

Polishers for fast removal
and contouring



Smooth - Choose the Red

Polishers for smoothing
and silky shine



Shine - Choose the Gray

Polishers for ultra
high shine

Order #	Kit	Quantity	Price
FPC01	Composite Polisher Complete Kit: 60 count	60 (10 each)	\$139.00
Order #	Refills		
FPCC010	Blue coarse grit point composite polisher	30	\$53.00
FPCC040	Blue coarse grit cup composite polisher	30	\$53.00
FPCM010	Red medium grit point composite polisher	30	\$53.00
FPCM040	Red medium grit cup composite polisher	30	\$53.00
FPCF010	Gray fine grit point composite polisher	15	\$53.00
FPCF040	Gray fine grit cup composite polisher	15	\$53.00



Do not exceed
maximum
speed



Wear
breathing
protection



Wear eye
protection



Use waterspray
minimum 50ml/
min



**INTRODUCTORY
PRICE
\$139**

For more information go to www.garrisdental.com. Be sure to check out the website to see
the full line of innovative Garrison products and our 6-month money-back guarantee.



Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

Offer expires 12/31/2017
To order: Call toll free **888.437.0032**
or online at www.garrisdental.com



Mention Code
IN1116RP

© 2016 Garrison Dental Solutions, LLC