

# Fit Strip™

**F**lexible  
**I**nterproximal finishers with a  
**T**wist



**FitStrips** are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

**One handy handle!** FitStrips' unique attachable handle has been **IMPROVED!** It's even easier to attach (or remove) and makes a world of difference for both clinician ease-of-use and patient comfort. Your hand stays outside the mouth, providing a clear field of view and easier operation.

**Do the Twist!** Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place—curved for finishing or straight for IPR. It doesn't get any easier than that!



**FitStrip™ makes your procedures easier!**

- Orthodontics: The first truly complete IPR/tooth slenderizing kit
- Composite restorations: contouring/finishing all interproximal surfaces
- Crown and bridge: cement removal and clean-up

**Ready for a marathon!** FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.



**F**lexible  
**I**nterproximal finishers with a  
**T**wist



## Your Complete IPR Solution!

**NEW!** .40mm and .50mm strips easily accommodate maximum reduction cases

**NEW!** IPR feeler gauges for exact reduction

**IMPROVED!** Attachable handle for greater ease of use and patient comfort

FitStrip™

**Kits** All kits come in a package of 10 FitStrips with 2 handles.

### NEW Universal Kit

1 IPR gauge  
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)  
6 double-sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK04 \$185

### NEW IPR Kit

1 IPR gauge  
10 double-sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

FPSK05 \$185

### Starter Kit

2 serrated  
4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)  
4 double-sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)

FPSK01 \$169

### Single Sided Kit

2 serrated  
8 single-sided (2 of each .08 mm, .10 mm, .13 mm, .18 mm)

FPSK02 \$169

### Double Sided Kit

2 serrated  
8 double-sided (2 of each .11 mm, .15 mm, .21 mm, .30 mm)

FPSK03 \$169

**NEW UNIVERSAL AND IPR  
FITSTRIP KIT  
\$185**






### IPR Kit

FPSK05 contains:

2 handles  
1 IPR gauge  
10 Strips (2 each: .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

**NEW!**



Refills		Single Sided	Double Sided	10 pk w/2 handles	4 ct packs (add – 4 to model #)
Super fine 15 Micron		FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick	Either Refill \$169	\$81
Fine 30 Micron		FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick	Either Refill \$169	\$81
Medium 46 Micron		FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick	Either Refill \$169	\$81
Coarse 90 Micron		FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick	Either Refill \$169	\$81
Medium 46 Micron			FPSTMDS 0.40 mm thick	Refill \$169	\$81
Medium 46 Micron			FPSXTMDS 0.50 mm thick	Refill \$169	\$81
Serrated		FPSSAW 0.05 mm thick		\$169	\$81
FitStrip replacement handles (2-pack) FPSHANDLE-2				\$29	
Gauges FPSGAUGE 				\$39	

For more information go to [www.garrisondental.com](http://www.garrisondental.com). Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.



**Garrison**  
Dental Solutions

150 DeWitt Lane  
Spring Lake, MI USA  
49456

Offer expires 12/31/2019  
To order: Call toll free **888.437.0032**  
or online at [www.garrisondental.com](http://www.garrisondental.com)

Mention Code  
**IN119FS**