


NEW!

NiTⁱⁿ

NiTⁱⁿ Sectional Matrix System


The new cost-effective solution for tight contacts from Garrison!

As the name implies, the new NiTin™ sectional matrix system separator rings are made from a special formulation of NiTinol nickel titanium for super strength and resiliency. Rather than typical material used by other brands, NiTin™ is made with drawn-wire Nitinol which produces superior spring characteristics. Add in the PEEK® reinforcement on the bow of the ring and you've got a system that produces exceptional contacts for hundreds and hundreds of uses.



Drawn-wire nickel titanium for the ultimate in long life!

Proprietary heat-set processing increases strength and resiliency.



PEEK® (polyether ether ketone) super-polymer reinforcement reduces ring stretch.

Ultra-durable PEEK® tips are the strongest, longest lasting polymer tips on any ring.

Precisely engineered tip geometry allows NiTin™ to grip where others pop-off—all while reducing flash.

Smaller, less "bulky" design fits well in tight spots and is still stackable for MODs.




All Garrison products are covered by our six-month money-back guarantee. If you are not satisfied with your product, return the unused portion and we will issue a refund. It's that simple.


Garrison
Dental Solutions



The Complete NiTin™ System



Full contour matrix bands for ideal anatomy!



Similar in dimension to the top-of-the-line Composi-Tight® 3D Fusion™ matrix bands, NiTin™ bands are slightly stiffer to improve ease of placement and are an economical choice for a full curve band.

50 ct packs \$60 each

NiTin's drawn-wire nickel titanium rings are strong and resilient!



NTG400 Long Rings (black) are excellent for tall molars and for easy stacking on MODs


\$195 pkg of 2 rings



NTG500 Standard Rings (white) are perfect for shorter teeth and next to rubber dam clamps.

\$195 pkg of 2 rings

Classic G-Wedges™ for a gentle yet firm gingival seal!



The G-Wedge is the perfect complement to the simple and cost effective NiTin™ system.

100 ct standard packs \$43 each

Precision Instrumentation completes the NiTin™ system!



\$75

Low profile tips simplify stacking on MODs with NiTin's smaller rings.



\$23


This surgical tissue forceps with its deep interdigitating teeth provide ample grip and tactile control for both band and wedge placement.



Starter Kit NTG-KMG-01

- 1 Long NiTin™ Ring
- 1 Standard NiTin™ Ring
- 100 Assorted NiTin™ Matrices
- 100 Assorted G-Wedges™
- 1 NiTin™ Ring Placement Forceps
- 1 NiTin™ Band Forceps

\$440



For more information go to www.garrisdental.com. Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.

Garrison
Dental Solutions

150 DeWitt Lane
Spring Lake, MI USA
49456

Offer expires 12/31/2019
To order: Call toll free **888.437.0032**
or online at www.garrisdental.com



Mention Code
IN719NT

© 2019 Garrison Dental Solutions, LLC