

BAND COMPARISON CHART

Garrison

SG Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Strata-G™ matrix bands with optimized marginal ridge						
thickness 0.0016"			NA			
height	SG050 3.5 mm	SG100 4.5 mm	NA	SG175 5.5 mm	SG200 6.5 mm	SG300 7.5 mm

Item	Standard (-M) Qty.	Bulk Qty.
SG050	50	100
SG100	50	100
SG175	50	100
SG200	50	100
SG300	50	100

FXH Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.						
thickness 0.0015"	NA					
height	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm

Item	Standard (-M) Qty.	Bulk Qty.
FXH100	50	100
FXH150	50	100
FXH175	50	100
FXH200	50	100
FXH300	50	100

FX Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.						
thickness 0.0016"	NA					
height	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm

Item	Standard (-M) Qty.	Bulk Qty.
FX100	50	100
FX150	30	60
FX175	50	100
FX200	50	100
FX300	30	60

SXR Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.						
thickness 0.0016"	NA					
height	NA	SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm w/ext 8.6 mm

Item	Standard (-M) Qty.	Price	Bulk Qty.
SXR100	50		100
SXR150	30		60
SXR175	50		100
SXR200	50		100
SXR300	30		60

SM Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Slick Bands™ Matrices have non-stick coating.						
thickness 0.0016"						
height	SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm


Item	Standard (-M) Qty.	Bulk Qty.
SM050	NA	60
SM100	50	100
SM150	30	60
SM175	50	100
SM200	50	100
SM300	30	60

M Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.						
thickness 0.0015"						
height	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm


Item	Standard (-M) Qty.	Bulk Qty.
M050	25	50
M100	50	100
M150	25	50
M175	50	100
M200	50	100
M300	25	50

AU Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
The Gold bands are the same as the Composi-Tight® bands.						
thickness 0.0015"				NA		
height	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm

Item	Standard (-M) Qty.	Bulk Qty.
AU050	25	50
AU100	50	100
AU150	25	50
AU200	50	100
AU300	25	50

B Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
The original bands have a slightly shorter length.						
thickness 0.0015"				NA		
height	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm

Item	Standard (-M) Qty.	Bulk Qty.
B050	25	50
B100	50	100
B150	25	50
B200	50	100
B300	25	50