










| | Prior Systems | All-New Composi-Tight® 3D Fusion™ |
|----------------|---|--|
| RINGS |  |  <p>All-New Separator Rings</p> <ul style="list-style-type: none"> • BOTH Short ring (blue) and Tall ring (orange) have Ultra-Grip™ retention extensions to prevent “spring-off” • Improved stacking for MODs • 20%+ Improved separating force • Enhanced marginal ridge built into Soft-Face™ silicone tips for better adaptation |
| WIDE PREP RING | <p>Not Available</p> |  <p>Wide Preparation Ring</p> <ul style="list-style-type: none"> • Dentistry's first Wide Preparation ring! Greatly simplifies these tougher cases – even with a missing cusp! • The same great retention and adaptation as the Short ring (blue) and Tall ring (orange) above |
| MATRIX BANDS |  |  <p>Full Curve Matrix Bands</p> <ul style="list-style-type: none"> • Recreates ideal anatomy with new Full Curve design • Easier placement with new Grab-Tabs™ on the occlusal edge • Longer bands are a perfect match for the Wide Prep ring |
| WEDGES |  |  <p>Ultra-Adaptive Wedges</p> <ul style="list-style-type: none"> • Soft-Face™ exterior adapts to interproximal irregularities • Flexible fins glide in but prevent wedge back-out • Superior gingival seal |
| FORCEPS |  |  <p>Dentistry's Premier Forceps</p> <ul style="list-style-type: none"> • Drop-forged durability prevents leverage-robbing flexion • Easily opens strong rings regardless of hand strength or size • Fits virtually any ring from any manufacturer |